



CRYSTAL AIRWAY MASK

- Soft and non-traumatic laryngeal protector avoid lesions and ensure a good hold and ventilation
- The rear cuff is reinforced with a thicker wall so that the device can be easily inserted
- Inflation line with valve is integrated into the wall of the ventilation tube
- **LATEX DEHP and BFA FREE**
- Available in 7 sizes
- Each size is color coded
- Sterile gas ETO
- Packaging Blister
- Single patient use

51.1-51.5 Distalinė dalis gerai prisitaiko ir priglunda prie viršutinių kvėpavimo takų, jos paviršius švelnus

51.1-51.5 Vožtuvas ir distalinės dalies išpūtimą rodantis balionas yra atokiau nuo vamzdelio konektoriaus

51.1-51.5 Distalinės dalies forma klasikinė Braino

51.1-51.5 Vamzdelis, skirtas distalinei daliai

51.1-51.5 Kvėpavimo vamzdelis permatomas

51.1-51.5 Vartotojui reikalinga informacija yra matomoje kvėpavimo vamzdelio dalyje



51 p.d.



• CRYSTAL AIRWAY MASK



kodas dydis

	Code	Size	Patient's weight Kg	I.D. internal	O.d. external	Maximum volume of the cuff	Items per box
	KM-877	1	< 5	5,6	8,6	5	10
51.1 p.d.	KM-867	1,5	5-10	6,4	10	7	10
51.2 p.d.	KM-857	2	10-20	7,8	12	10	10
	KM-847	2,5	20-30	9,5	14	15	10
51.3 p.d.	KM-837	3	30-50	10	16	20	10
51.4 p.d.	KM-827	4	50-70	11,4	18	35	10
51.5 p.d.	KM-817	5	➤ 70	11,9	20	48	10