

Katalog-Nr. Kategorie		Bezeichnung	MERLIN Diagnostika
E1-199-040 H		MICRONAUT-S Pneumococci & Haemophilus MIC	
Layout	K		
1 Test			

	1	2	3	4	5	6	7	8	9	10	11	12
A	COL 2	AMP 16	IMP 4	MER 8	CRO 16	CXM 16	ERY 16	CLR 32	CIP 8	RAM 0,25	RAM 2	DOX 16
B	PEN 4	AMP 8	IMP 2	MER 4	CRO 8	CXM 8	ERY 8	CLR 16	CIP 4	RAM 0,125	RAM 1	DOX 8
C	PEN 2	AMP 4	IMP 0,5	MER 2	CRO 4	CXM 4	ERY 4	CLR 8	CIP 2	RAM 0,0625	RAM 0,5	DOX 4
D	PEN 1	AMP 2	IMP 0,125	MER 1	CRO 2	CXM 2	ERY 2	CLR 4	CIP 1	MOX 2	CLI 2	DOX 2
E	PEN 0,5	AMP 1	AMC 4/2	MER 0,5	CRO 1	CXM 1	ERY 1	CLR 2	CIP 0,5	MOX 1	CLI 1	DOX 1
F	PEN 0,25	AMP 0,5	AMC 2/2	MER 0,25	CRO 0,5	CXM 0,5	ERY 0,5	CLR 1	CIP 0,25	MOX 0,5	CLI 0,5	DOX 0,5
G	PEN 0,125	AMP 0,25	AMC 1/2	MER 0,125	CRO 0,25	CXM 0,25	ERY 0,25	CLR 0,5	CIP 0,125	MOX 0,25	CLI 0,25	DOX 0,25
H	PEN 0,0625	AMP 0,125	AMC 0,5/2	MER 0,0625	CRO 0,125	CXM 0,125	ERY 0,125	CLR 0,25	CIP 0,0625	MOX 0,125	CLI 0,125	GC

[Angabe der Konzentrationen in µg/ml]

[Ausgabe von DIN-Kürzeln]

LEGENDE	
AMC	Amoxicillin/Clavulansäure
AMP	Ampicillin
CRO	Ceftriaxon
CXM	Cefuroxim
CIP	Ciprofloxacin
CLR	Clarithromycin
CLI	Clindamycin
COL	Colistin
DOX	Doxycyclin
ERY	Erythromycin
IMP	Imipenem
MER	Meropenem
MOX	Moxifloxacin
PEN	Penicillin G
RAM	Rifampicin
GC	Wachstumskontrolle