

16<sup>th</sup> October 2025

## STATEMENT LETTER

TO WHOM IT MAY CONCERN

We, SCW Medicath, hereby state that the Syringes 832624 / 42.16.40000 use with the MEDTRON Accutron HP injector:

- syringe - 200 ml;
- sterile (symbol on the package);
- disposable (marked with a symbol);
- the syringe mounting connection is adapted to the automatic MEDTRON Accutron HP injector;
- the pin is adapted for the injection of contrast into a syringe or a quick-filling tube;
- high-pressure transparent (transparent), soft, flexible, spiral extension line 150 cm, for connecting the syringe with the intravenous catheter;
- confirming that the proposed kit is suitable for work with MEDTRON Accutron HP injector;
- packed in one set;
- the expiration date of the product is marked on the package;
- with a designated place to open the package.

Yo

SCW MEDICATH LTD

SCW MEDICATH LTD  
Steve Wang (Sales Manager)

深圳市益心达医学新技术有限公司

Signature: