

# Knitted vascular prosthesis with collagen Ra K linear

CE  
1434

## KNITTED VASCULAR PROSTHESIS OFFERS:

### Direct use

The collagen layer allows the prosthesis to be used without the need for any modification before implantation.

### Easy handling

Coloured guiding line for easy orientation in the surgical wound and easy prosthesis handling during the surgery.

### Excellent features

The unique design ensures the strength of the prosthesis wall while providing the required porosity, pliability, softness, and ease of sewing.

### Simple adjustment

The prosthesis can be cut to the desired shape at the ends with scissors, the edges do not fray and thus allow the creation of a fine anastomosis.

### Parameters

Water permeability and leakage  $<5 \text{ ml.cm}^{-2}.\text{min}^{-1}$ . Impregnation with crosslinked bovine collagen, formaldehyde and glutaraldehyde free.

### Indications

The knitted vascular prosthesis with collagen is indicated for surgical treatment of aneurysms and arterial occlusions. The prosthesis is especially suitable for **the reconstruction of the abdominal aorta and peripheral arteries.**

## ORDERING INFORMATION (LINEAR)

Ø [mm]	Catalogue number length 100 mm*	Catalogue number length 200 mm*	Catalogue number length 250 mm*	Catalogue number length 300 mm
6	022-06/100	022-06/200	022-0061	022-06/300
7	022-07/100	022-07/200	022-0071	022-07/300
8	022-08/100	022-08/200	022-0081	022-08/300
9	022-09/100	022-09/200	022-0091	022-09/300
10	022-10/100	022-10/200	022-0101	022-10/300
12	022-12/100	022-12/200	022-0121	022-12/300
14	022-14/100	022-14/200	022-0141	022-14/300
16	022-16/100	022-16/200	022-0161	022-16/300
18	022-18/100	022-18/200	022-0181	022-18/300
20	022-20/100	022-20/200	022-0201	022-20/300
22	022-22/100	022-22/200	022-0221	022-22/300
24	022-24/100	022-24/200	022-0241	022-24/300
26	022-26/100	022-26/200	022-0261	022-26/300
28	022-28/100	022-28/200	022-0281	022-28/300
30	022-30/100	022-30/200	022-0301	022-30/300
32	022-32/100	022-32/200	022-0321	022-32/300
34	022-34/100	022-34/200	022-0341	022-34/300

\* These items are not a standard item, produced on the customer's request.

Ø [mm]	Catalogue number length 400 mm*	Catalogue number length 500 mm	Catalogue number length 600 mm*	Catalogue number length 700 mm*
6	022-06/400	022-0062	022-0063	022-0064
7	022-07/400	022-0072	022-0073	022-0074
8	022-08/400	022-0082	022-0083	022-0084
9	022-09/400	022-0092	022-0093	022-0094
10	022-10/400	022-0102	022-0103	022-0104
12	022-12/400	022-0122	022-0123	022-0124
14	022-14/400	022-0142	022-0143	022-0144
16	022-16/400	022-0162	022-0163	022-0164
18	022-18/400	022-0182	022-0183	022-0184
20	022-20/400	022-0202	022-0203	022-0204
22	022-22/400	022-0222	022-0223	022-0224
24	022-24/400	022-0242	022-0243	022-0244
26	022-26/400	022-0262	022-0263	022-0264
28	022-28/400	022-0282	022-0283	022-0284
30	022-30/400	022-0302	022-0303	022-0304
32	022-32/400	022-0322	022-0323	022-0324
34	022-34/400	022-0342	022-0343	022-0344

\* These items are not a standard item, produced on the customer's request.



Sterilized by irradiation  
Expiration 3 years from the date of manufacture