

Debridement

Wound bed preparation is a key element to ensuring optimum wound healing. Wound debridement removes barriers such as non-viable tissue, biofilm and foreign material and/or debris from the wound bed.

Mechanical is the most common form of debridement as it is quick, safe and effective with no specialist training required.

Kliniderm® debride is designed to provide rapid, gentle, mechanical debridement. Its intuitive design is double sided with a continuous soft edge, the product is placed over a gloved hand or finger enabling improved handleability and control during debridement.

Klinion® Kliniderm debride

Monofilament polyester mechanical debrident. Removes exudate, cellular debris and keratosis. Gentle on intact tissue. Medical Device. Sterile if pouch and seal are intact.

ART	REF	Size	Pieces inner carton	Pieces outer carton
3059040	40511800	14 x 10 CM	5	100
3059041	40511801	8 x 4 CM	5	60



Skin adhesive

Klinion® skin adhesive

Topical skin adhesive (SurgiSeal Stylus®). Professional liquid skin adhesive containing a monomeric (2-octyl cyanoacrylate) formulation and the colourant D & C Violet #2. Klinion® skin adhesive is intended for topical applications only to hold closed easily approximated skin edges of wounds from surgical incisions including punctures from minimally invasive surgery and simple, thoroughly cleansed, trauma induced lacerations. Klinion® skin adhesive may be used in conjunction with, but not in place of, deep dermal sutures. Medical Device. Sterile if pouch and seal are intact.

ART	REF	Size	Pieces inner carton	Pieces outer carton
3063854	40115500	0.5 ML	12	576



Additional advanced dressings

A range of advanced dressings containing alginate, silicone wound contact layer, hydrocolloid, films and Klintulle (paraffine dressing), as well as a range of honey products (Mestran).

Klinion® Kliniderm hydrocolloid

Hydrocolloid dressing for superficial wounds with low exudate, does not adhere to the wound, self-adhesive. Sterile.

ART	REF	Size	Pieces inner carton	Pieces outer carton
3036409	174302	10 x 10 CM	5	100
3036410	174303	15 x 15 CM	5	50
3036412	174304	20 x 20 CM	5	50



Klinion® Kliniderm hydrocolloid border

Hydrocolloid dressing for superficial wounds with low exudate, does not adhere to the wound. The adhesive border enables an optimal fixation. Sterile.

ART	REF	Size	Pieces inner carton	Pieces outer carton
3036391	174201	7 x 7 CM	5	100
3036394	174202	10 x 10 CM	5	100
3036396	174203	14 x 14 CM	5	50



Klinion® Kliniderm hydrocolloid border sacrum

Hydrocolloid dressing for superficial wounds with low exudate, does not adhere to the wound. The adhesive border enables an optimal fixation. Sterile.

ART	REF	Size	Pieces inner carton	Pieces outer carton
3036398	174204	15 x 18 CM	5	50



Klinion® Kliniderm hydrocolloid thin

Hydrocolloid dressing for superficial wounds with low exudate, does not adhere to the wound, self-adhesive. Sterile.

ART	REF	Size	Pieces inner carton	Pieces outer carton
3036414	174101	7.5 x 7.5 CM	5	300
3036357	174102	10 x 10 CM	5	300
3036358	174103	15 x 15 CM	5	150
3036359	174104	20 x 20 CM	5	100

