



## DOCON CE1275

State-of-the-art blood collection mixer

**DESIGNED  
TO EASE  
YOUR WORK**



### BENEFITS

- Large bag tray, smart tube detector, large display, remote control with adjustable height.
- Long life batteries allowing up to 1 week or 70 donations without recharge.
- High weighing precision (< 1%) and efficient mixing.
- High level of traceability through a large variety of data exchange features.
- Fully customizable collection programs.
- Integrated tube sealer (MacoSeal) as an option.



### SEALING FUNCTION/MACOSEAL

Maintains the integrity of the donation. The Docon is equipped with the MACOSEAL quality hand sealer. Sealing operations can be recorded in the Docon with the donation data.



Docon<sup>®</sup>, *Donation Controller*, is a **FULLY AUTOMATIC** blood collection mixer that guarantees **OPTIMAL** blood collection and maintains the requirements for **HIGH QUALITY** blood components.

Large bag tray to receive all types of blood bags.

Very large display (9 cm x 12 cm) : information on the current status of the donation (duration of the donation, date, time and battery status). Measurement of blood flow and a flow graphic can be displayed upon demand.

**Sensoclamp**  
Equipped with a smart tube recognition to guarantee the highest security during collection. Automatic clamping upon insertion of the tube.

Robust transport case.



**Docon Mobile II**  
Reduces the need to bend down. Shows the current status of the donation. Provided with a keyboard : manual input of information such as temperature, haemoglobin etc. Can be operated from a flexible gooseneck holder.

- Very quick start-up
- Fully customized and adapted to each specific need. Ex: type of bags, different sounds and visual alarms, bar code controls.

Go to [www.macopharma.com](http://www.macopharma.com) to get a complete overview of our products.

DOCON is CE certified, Class IIa MDD, 93/42/EEC.

MACO PHARMA  
Rue Lorthiois  
F-59420 Mouvaux (France)  
[www.macopharma.com](http://www.macopharma.com)

Distributed by:

