

# Fistula needles

## Features and benefits

**Smooth and painless puncture**  
Specially cut bevel design  
New dry siliconization method  
Ultrathin 95 µm wall

### Optimal blood flow characteristics

Sidehole on arterial needles  
Specially-designed Y connector

### Simplified handling

Rotatable colour-coded wings  
Colour-coded indication of the cannula position  
Luer connectors for easy and safe connection  
Soft anti-kinking PVC tubing and easy clamps


### Packaging

Inner box : 50 needles  
Outer box : 500 needles

### Gamma sterilized

### Gambro Plaster XL:

cellulose sponge expands rapidly in contact with fluids, resulting in immediate pressure on the wound. Can be removed easily when hemostasis has been achieved.



The size of the fistula needle must always suit the chosen blood flow.

**GAMBRO®**





	Cannula			Tubing length mm	Arterial with sidehole	Venous w/o sidehole	Wings	
	Gauge	diameter mm	length mm				fixed	color
8501410000	14	2,0	25	150		○		lavender
8501410100	14	2,0	25	150	○			lavender
8501420000	14	2,0	25	300		○		lavender
8501420100	14	2,0	25	300	○			lavender
8501510000	15	1,8	25	150		○		brown
8501510100	15	1,8	25	150	○			brown
8501511100	15	1,8	25	150	○		X	brown
8501520000	15	1,8	25	300		○		brown
8501520100	15	1,8	25	300	○			brown
8501521000	15	1,8	25	300		○	X	brown
8501610000	16	1,6	25	150		○		green
8501610100	16	1,6	25	150	○			green
8501611100	16	1,6	25	150	○		X	green
8501620000	16	1,6	25	300		○		green
8501620100	16	1,6	25	300	○			green
8501621000	16	1,6	25	300		○	X	green
8501621100	16	1,6	25	300	○		X	green
8501710000	17	1,5	25	150		○		red
8501710100	17	1,5	25	150	○			red
8501720000	17	1,5	25	300		○		red
8501720100	17	1,5	25	300	○			red
8511410100	14	2,0	25	150	Single-Needle			lavender
8511510100	15	1,8	25	150	Single-Needle			brown
8511520100	15	1,8	25	300	Single-Needle			brown
8511610100	16	1,6	25	150	Single-Needle			green
8511620100	16	1,6	25	300	Single-Needle			green
8521610100	16	1,6	25	150	○	○		green
8N02057002	Plaster XL	Bandage: 30x72 mm			Sponge: 15x20 mm			

Other models (e.g. with different cannula length) might be available on request.

For further information and operating instructions, please refer to operator's manual.

Fistula needle — FSN 15 AS — S=Sidehole

SN=Single needle — Tube length:

Cannula diameter in Gauge — A=150 mm B=300 mm

**Fistulės adatos**  
**Ypatybės ir pranašumai**

16G  
15G  
14G

100    200    300    400    500  
kraujo srovė (ml/min.)

Fistulės adatos dydis visada turi atitikti pasirinktą kraujo srovės greitį.

**GAMBRO**

*/Vertimas iš anglų kalbos/*

**Sklandus ir neskausmingas dūris**

Speciali įstrižo pjūvio forma  
Naujas dengimo sausu silikonu  
metodas

Ypač plona 95 µm sienelė

**Tinkamiausios kraujo srovės  
charakteristikos**

Šoninė anga arterijoms skirtose  
adatosė

Specialiai sukurta „Y“ jungtis

**Paprastesnis naudojimas**

Spalvomis koduojami pasukamieji  
sparneliai

Spalvomis koduojama kaniulės  
padėties nuoroda

„Luer“ tipo jungtys, užtikrinančios  
nesudėtingą ir saugų prijungimą

Minkšti ir neužsilenkiantys PVC  
vamzdeliai ir lengvai valdomi  
spaustukai

**Pakuotė**

Vidinė dėžė: 50 adatų

Išorinė dėžė: 500 adatų

**Sterilizuota gama spinduliais**

**„Gambro Plaster XL“:**

celiuliozės kempinė, kuri susilietusi  
su skysčiais greitai išsiplečia ir  
prispaudžia žaizdą. Sustabdžius  
kraujavimą ją galima lengvai nuimti.



		Kaniulė			Vamzdelių ilgis mm	Arterinis su šonine anga	Veninis be šoninės angos	Sparneliai	
		skersmuo		ilgis mm				fiksuota	spalva
		dydis	mm						
8501410000	F 14A	14	2,0	25	150		o		levandų
8501410100	F 14AS	14	2,0	25	150	o			levandų
8501420000	F 14B	14	2,0	25	300		o		levandų
8501420100	F 14BS	14	2,0	25	300	o			levandų
8501510000	F 15A	15	1,8	25	150		o		ruda
8501510100	F 15AS	15	1,8	25	150	o			ruda
8501511100	F	15	1,8	25	150	o		X	ruda
8501520000	F 15B	15	1,8	25	300		o		ruda
8501520100	F 15BS	15	1,8	25	300	o			ruda
8501521000	F 15BF	15	1,8	25	300		o	X	ruda
8501610000	F 16A	16	1,6	25	150		o		žalia
8501610100	F 16AS	16	1,6	25	150	o			žalia
8501611100	F	16	1,6	25	150	o		X	žalia
8501620000	F 16B	16	1,6	25	300		o		žalia
8501620100	F 16BS	16	1,6	25	300	o			žalia
8501621000	F 16BF	16	1,6	25	300		o	X	žalia
8501621100	F	16	1,6	25	300	o		X	žalia
8501710000	F 17A	17	1,5	25	150		o		raudona
8501710100	F 17AS	17	1,5	25	150	o			raudona
8501720000	F 17B	17	1,5	25	300		o		raudona
8501720100	F 17AS	17	1,5	25	300	o			raudona
8511410100	FSN	14	2,0	25	150	Viena adata			levandų
8511510100	FSN	15	1,8	25	150	Viena adata			ruda
8511520100	FSN	15	1,8	25	300	Viena adata			ruda
8511610100	FSN	16	1,6	25	150	Viena adata			žalia
8511620100	FSN	16	1,6	25	300	Viena adata			žalia
8521610100	F 16AS/A	16	1,6	25	150	o	o		žalia
8N02057002	„Plaster XL“	tvarstis: 30 x 72 mm				kempinė: 15 x 20 mm			

Paprašius galima įsigyti kitokių modelių (pvz., su kitokio ilgio kaniule).  
Daugiau informacijos ir naudojimo nurodymus rasite vartotojo vadove.

Fistulės adata FSN 15 AS, S – šoninė anga,  
SN – viena adata  
Vamzdelio ilgis:  
Kaniulės skersmuo: A – 150 mm, B – 300 mm

HCEN1390 © „Gambro Lundia AB“. 2002 m. gegužės mėn. Spausdino „Elanders Skogs Grafiska AB“,  
1 peržiūra, 2002 m. gegužės mėn.

Aš, vertėja Ana Bulienė,  
esu susipažinusi su LR BK 235 straipsniu,  
kuris numato baudžiamąją atsakomybę  
už neteisingą vertimą.  
Bulienė .....

