



STAINS & INSTRUMENTS FOR INTELLIGENT MICROSCOPY


KIT RAL 555

Fast Versatile Staining

HAEMATOLOGY




PARASITIOLOGY



HISTOLOGY



CYTO-BACTERIOLOGY



Kit RAL 555

Fast version of May-Grunwald Giemsa Staining

EASY TO USE

FAST

The RAL 555 Kit is very easy to use since it consists in successive baths with direct immersion of the slide in the vials.

- The color of the containers results in an enhanced visibility and in a more rapid identification of the various solutions.
- The vials, with broader necks, make it easier to insert the slide.

The kit packaging must not be discarded as it serves as a test stand.

EFFICIENT

- Sensitivity, reading and interpretation of the staining are similar to those of traditional M.G.G stainings.
- Ideal for emergency diagnosis thanks to its quick execution.
- Formulation and packaging have been optimised for the user confort and safety.



Field of Application	Processing time
Haematology :	
Blood smears	= 15 secondes
Bone marrow	= 3 minutes
Cyto-bacteriology :	
Liquid Cytologyc investigation (Needle punctures, urine)	= 15 secondes
Parasitology :	
Plasmodium detection test	= 1minute and 4 secondes

For ordering, please contact your distributor.

Also available in separate bottles : Kit RAL 555 : 361550-0000

FIX-RAL 555 : 362870-0100, -1000, -2500 ml

EOSINE-RAL 555 : 361640-0100, -1000, -2500 ml

BLEU-RAL 555 : 361550-0100, -1000, -2500 ml

RAL DIAGNOSTICS - Site Montesquieu - 33651 Martillac FRANCE

Tel : + 33 (0) 5 57 96 04 04 - Fax : +33 (0) 5 57 96 04 05

www.ral-diagnostics.com