



Vilnius

**Viešoji įstaiga Vilniaus universiteto ligoninė Santaros klinikos** (toliau – **Pirkėjas**), atstovaujamos generalinio direktoriaus Felikso Jankevičiaus, veikiančio pagal įstaigos įstatus ir **UAB Medfarmos laboratorijos** (toliau – **Pardavėjas**), atstovaujama direktorės Živilės Sukackienės veikiančios pagal įmonės įstatus, laimėjusi 2020 m. vasario 29 d. skelbtą atviro konkurso būdu vykusį „Diagnostinių reagentų, skirtų molekuliniais citogenetiniams tyrimams VULSK Nr. 20836“ viešąjį pirkimą (pirkimo skelbimo Nr. 2020-640919 ) (toliau - Pirkimas), sudarė šią sutartį:

## 1. SUTARTIES OBJEKTAS

1.1. Pardavėjas įsipareigoja parduoti Pirkėjui reagentus ir priemones skirtas molekuliniais citogenetiniams tyrimams (toliau – Prekės), o Pirkėjas įsipareigoja priimti jas ir už jas sumokėti šioje sutartyje nurodytą kainą sutarties 6 p. numatyta tvarka.

## 2. SUTARTIES VERTĖ

2.1. Sutarčiai taikomas **fiksuotų įkainių maksimalios sumos kainodaros metodas**. Sutartyje nustatyti **fiksuoti prekių įkainiai ir perkamų prekių preliminarūs kiekiai**, kurie nurodyti šios sutarties priede Nr. 1.

2.2. Pradinė (maksimali) sutarties vertė be PVM **33 333,34** Eur (trisdešimt trys tūkstančiai trys šimtai trisdešimt trys eurai 34 ct).

PVM suma 1 666,66 Eur (vienas tūkstantis šeši šimtai šešiasdešimt šeši eurai 66 ct).

**Pradinė (maksimali) sutarties vertė su PVM 35 000,00** Eur (trisdešimt penki tūkstančiai eurų 00 ct).

2.3. Įkainiai sutarties vykdymo laikotarpiu negalės būti keičiami per visą sutarties vykdymo laikotarpį, išskyrus kai pasikeičia pridėtinės vertės mokesčio (PVM). Perskaičiavimas vykdomas po Lietuvos Respublikos pridėtinės vertės mokesčio įstatymo, kuriuo keičiasi mokesčio tarifas, įsigaliojimo dienos. Pasikeitus PVM tarifo dydžiui, nepateiktų prekių įkainiai keičiami (mažinama ar didinama) proporcingai PVM pasikeitusio tarifo dydžiui. Įkainių pakeitimas įforminamas papildomu rašytiniu šalių susitarimu.

2.4. Į Prekių įkainius įskaityti visi Pardavėjo mokami mokesčiai, prekės pristatymo ir kitos išlaidos.

2.5. **Sutarties kaina** – už prekes pagal Sutartį tiekėjo gaunama ekonominė nauda.

2.6. Pirkimo sutarties sąlygos sutarties galiojimo laikotarpiu negali būti keičiamos, išskyrus tokias pirkimo sutarties sąlygas, kurias pakeitus nebūtų pažeisti Viešųjų pirkimų įstatymo 17 straipsnyje nustatyti principai ir tikslai, bei esant Viešųjų pirkimų įstatymo 89 straipsnyje nustatytioms sąlygoms.

2.7. Pirkėjas įsigys techninėje specifikacijoje išvardintus diagnostikos reagentus ir zondus, skirtus formalinu fiksuotiems ir parafinu impregnuotiems (toliau – FFPI) mėginiams arba ląstelių preparatams, in situ hibridizacijos (toliau – ISH) tyrimams atlikti. Išperkamų prekių kiekis gali kisti, priklausomai nuo molekulinų citogenetinių tyrimų poreikio: įsigijamų prekių kiekis  $\pm 30$  proc. nuo konkurso laimėtojo pasiūlyme kiekvienam tyrimui/reagentui siūlomo preliminarus prekių kiekio, tačiau neviršijant pirkimui numatytos pradinės (maksimalios) sutarties vertės be PVM.

2.8. Pirkėjas, atskirai kiekvienai pirkimo daliai, gali įsigyti tiekėjo sąrašė nenurodytų, tačiau su pirkimo objektu susijusių prekių (pavyzdžiui, naujesnė to paties reagento versija) neviršijant 10 procentų pradinės (maksimalios) sutarties vertės. Už prekių sąrašė nenurodytas, tačiau su pirkimo objektu susijusias prekes ir (ar) paslaugas bus apmokėta ne didesnėmis nei užsakymo dieną tiekėjo prekybos vietoje, kataloge ar interneto svetainėje nurodytomis galiojančiomis šių prekių ir (ar) paslaugų kainomis arba, jei tokios kainos neskelbiamos, tiekėjo pasiūlytomis, konkurencingomis ir rinką atitinkančiomis kainomis.

### 3. PERKAMOS PREKĖS

3.1. Šia Sutartimi Pardavėjas įsipareigoja perleisti prekes, kurių detali specifikacija ir įkainiai nustatyti šios Sutarties priede Nr.1. **Sutarties vykdymo metu įsigijami prekių kiekiai priklauso nuo faktinių Pirkėjo užsakymų, tačiau negali būti viršyta šios sutarties 2.2. punkte nurodyta maksimali suma, išpirkti mažesni kiekį Pirkėjas gali.**

3.2. Jei dėl nuo Pardavėjo nepriklausančių aplinkybių, kurių nebuvo įmanoma numatyti rengiant pirkimo dokumentus ir/ar sutarties sudarymo metu, Pardavėjas negali pristatyti pasiūlyme nurodyto modelio prekių, sutarties šalims išreiškus sutikimą, nekeičiant sutarties kainos, Pardavėjas gali pristatyti kito modelio prekes su sąlyga, kad naujas modelis visiškai atitiks pirkimo dokumentuose ir sutarties priede keliamus reikalavimus ir bus pristatomas už tą pačią kainą.

### 4. PREKIŲ PRISTATYMO TVARKA

4.1. Pardavėjas įsipareigoja pristatyti prekes savo transportu ne vėliau kaip per 14 (keturiolika) darbo dienų nuo Pirkėjo užsakyme nurodytos datos, o esant ypatingiems atvejams - per 7 (septynias) darbo dienas.

4.2. Prekių pristatymo vieta – Valstybinis patologijos centras, viešosios įstaigos Vilniaus universiteto ligoninės Santaros klinikų filialas (P. Baublio g. 5, LT-08406 Vilnius).

4.3. Pardavėjo paskirtas atsakingas asmuo: Živilė Sukackienė, tel.: 8 5 270 01 01, [info@medfarmos.lt](mailto:info@medfarmos.lt);

4.4. Pirkėjo paskirtas atsakingas asmuo: Aida Laurinavičienė, tel. +370 5247 7359 el. paštas: [Aida.Laurinaviciene@vpc.lt](mailto:Aida.Laurinaviciene@vpc.lt));

4.5. Paskirtas atsakingas asmuo už sutarties ir pakeitimų paskelbimą – viešųjų pirkimų specialistas (Algimantas Varžgalys, tel.+370 52477365, el. paštas: [algimantas.varzgalys@vpc.lt](mailto:algimantas.varzgalys@vpc.lt)).

### 5. PREKIŲ KOKYBĖ IR GARANTIJA

5.1. Pardavėjas garantuoja, kad parduodamos Prekės yra be defektų ir jų kokybė, žymėjimas, informacija vartotojui atitinka Pirkimo dokumentų nuostatas, galiojančius standartus, technines sąlygas ar kitus norminius aktus.

5.2. Pardavėjas pristato prekes tvarkingose pakuotėse. Pristatant prekes Pardavėjas pateikia prekės vartotojo instrukcijas lietuvių kalba (arba/ir ..... kalba, jei tai nustatyta pirkimo sąlygose). Prekių žymėjimas ant pakuotės turi būti lietuvių kalba (jei prekės gamintojo nėra žymimos valstybine kalba – pasitelkiant lipdukus ar kt. priemones).

5.3. Pateiktų prekių galiojimo terminas pristatymo metu turi būti ne trumpesnis nei 12 mėn. nuo gamintojo nustatyto galiojimo termino pabaigos, išskyrus atvejus, kai gamintojas nustato trumpesnį prekių galiojimo terminą.

5.4. Jei Pardavėjas Pirkėjui pateikia nekokybiškas, neatitinkančias šios sutarties sąlygų prekes, Pardavėjas privalo per 5 (penkias) darbo dienas nuo Pirkėjo pranešimo jas pakeisti kokybiškomis, atitinkančiomis šios sutarties sąlygas prekėmis.

5.5. Prekių realizacijos termino galiojimo metu pastebėjus, kad prekių kokybė neatitinka šios sutarties reikalavimų, po Pirkėjo pranešimo, ne ilgiau kaip per 5 (penkias) darbo dienas, Pardavėjas privalo pakeisti prekes.

### 6. ATSISKAITYMAS TARP ŠALIŲ

6.1. Pirkėjas apmoka Pardavėjui už prekes pagal gautas PVM sąskaitas faktūras per 30 dienų nuo prekių gavimo ir sąskaitos faktūros pateikimo dienos.

Jei mokėjimai pagal sutartis visiškai arba iš dalies atliekami iš tarpinių finansuojančių organizacijų gautomis lėšomis, taip pat kitomis objektyviai pagrįstomis aplinkybėmis, atsiskaitymo terminas Pirkėjo gali būti pratęstas iki 60 dienų nuo prekių gavimo ir sąskaitos faktūros pateikimo dienos.

6.2. Pridėtinės vertės mokesčio sąskaitos faktūros, sąskaitos faktūros, kreditiniai ir debetiniai dokumentai turi būti teikiami naudojantis informacinės sistemos „E. sąskaita“ priemonėmis. Mokėjimo dokumentų nepateikus „E. sąskaita“ priemonėmis, Pirkėjas turi teisę neatlikti mokėjimo.

## 7. SUTARTIES ĮVYKDYMO UŽTIKRINIMAS, ATSAKOMYBĖ

7.1. Sutarties įvykdymo užtikrinimo priemonė yra netesybos. Jei tiekėjas dėl savo kaltės vėluoja pristatyti Prekes, Pirkėjas turi teisę be rašytinio įspėjimo ir nesumažindamas kitų savo teisių gynimo priemonių, numatytų sutartyje, pradėti skaičiuoti delspinigius už kiekvieną vėluojamą prekių pristatymo dieną. Pardavėjo vėluojamų pristatyti prekių kaina mažinama 0,02 % nuo vėluojamų pristatyti prekių vertės už kiekvieną termino praleidimo dieną. Delspinigių suma gali būti išskaičiuojama iš Pardavėjui mokėtinų sumų.

7.2. Jei Pardavėjas iš esmės pažeidžia Sutartį 9.2.p. nurodytais atvejais, Pirkėjas turi teisę taikyti 10 % sutarties pradinės vertės baudą, įskaitant delspinigius.

7.3. Jei Pirkėjas dėl savo kaltės vėluoja atsiskaityti su Pardavėju per sutarties 6.1.p. numatytą terminą, jis įsipareigoja sumokėti 0,02 % dydžio delspinigius nuo nesumokėtos sumos.

## 8. FORCE MAJEURE

8.1. Force Majeure sąlygos taikomos vadovaujantis LR Vyriausybės 1996 m. liepos 15 d. nutarimu Nr. 840 patvirtintomis „Atleidimo nuo atsakomybės dėl nenugalimos jėgos (Force majeure) aplinkybėmis“, taisyklėmis.

## 9. SUTARTIES GALIOJIMAS IR NUTRAUKIMAS

9.1. Sutartis įsigalioja šalims ją pasirašius ir galioja iki visiško šalių įsipareigojimų pagal šią sutartį įvykdymo, bet ne ilgiau kaip **36 (trisdešimt šešis) mėnesius**.

9.2. Pirkėjas turi teisę vienašališkai nutraukti Sutartį įspėjęs Pardavėją prieš 15 darbo dienų, jeigu Pardavėjas ją iš esmės pažeidė:

9.2.1. parduota prekė yra netinkamos kokybės ir jos trūkumų neįmanoma pašalinti per protingą ir Pirkėjui priimtina terminą;

9.2.2. Pardavėjas prekių nepristatė per Sutartyje nustatytą bei Pirkėjo nurodytą papildomą protingą terminą.

9.3. Pardavėjas turi teisę vienašališkai nutraukti Sutartį įspėjęs Pirkėją prieš 15 darbo dienų, jeigu Pirkėjas ją iš esmės pažeidė:

9.3.1. Pirkėjas daugiau kaip du kartus laiku nesumokėjo už Prekes, kai jos buvo perduotos nustatytais terminais;

9.3.2. Pirkėjas daugiau kaip du kartus nepriėmė tinkamoms kokybės Prekių, kai jos buvo perduotos nustatytais terminais.

9.4. Sutartis taip pat gali būti nutraukta Šalių raštišku susitarimu.

## 10. KITOS SĄLYGOS

10.1. Pardavėjas įsipareigoja išrašomoje sąskaitoje-faktūroje vartoti tuos pačius prekių pavadinimus ir mato vnt., kokie yra pridedamoje specifikacijoje. Taip pat ant sąskaitos faktūros privalo užrašyti sutarties numerį ir datą, pagal kurią parduodamos prekės.

10.2. Sutarties vykdymo metu Pardavėjo gauta informacija ir dokumentai yra konfidencialūs. Be išankstinio raštiško Pirkėjo leidimo Pardavėjas neskelbia ir neatskleidžia jokių sutarties nuostatų, išskyrus atvejus, kai tai būtina vykdant sutartį.

10.3. Ginčai, kylantys tarp šalių, sprendžiami šalių derybomis, o nepavykus jų išspręsti – teismine tvarka Lietuvos Respublikos teismuose.

10.4. Sutartis sudaryta dviem egzemplioriais, po vieną šalims, lietuvių kalba. Abu egzemplioriai turi vienodą juridinę galią.

10.5. Pardavėjas patvirtina, kad į parduodamas prekes tretieji asmenys neturi jokių teisių.

10.6. Sudarius Sutartį, tačiau ne vėliau negu Sutartis pradėdama vykdyti, Pardavėjas įsipareigoja pranešti tuo metu žinomų subtiekiųjų pavadinimus, kontaktinius duomenis ir jų atstovus. Pirkėjas taip pat reikalauja, kad Pardavėjas informuotų apie minėtos informacijos pasikeitimus visu Sutarties vykdymo metu, taip pat apie naujus subtiekiuosius, kuriuos jis ketina pasitelkti vėliau. /Pasiūlyme Pardavėjas nenurodė, kad pasitelks subtiekiuosius/.

10.7. Pardavėjas gali keisti Sutartyje nurodytus subtiekiuosius tik prieš tai raštu pranešęs Pirkėjui apie tokio keitimo būtinybę ir gavęs jo raštišką sutikimą.

10.8. Pardavėjas Sutarties vykdymo metu gali inicijuoti subtiekiėjo, numatyto Sutartyje, pakeitimą, nurodydamas tokio keitimo motyvus.

10.9. Pirkėjui sutikus su subtiekiėjo pakeitimu, Pirkėjas kartu su Pardavėju raštu sudaro susitarimą dėl subtiekiėjo pakeitimo, kurį pasirašo Šalys. Šis susitarimas yra neatskiriama Sutarties dalis PRIDEDAMA. Specifikacija ir įkainiai.

## 1. ŠALIŲ ADRESAI IR REKVIZITAI

### Pardavėjas

UAB „Medfarmos laboratorijos“  
P. Smuglevičiaus g. 1, Vilnius LT-08311  
Įmonės kodas 124484379  
PVM kodas LT244843716  
Tel. 8 5 270 01 01, el. paštas info@medfarmos.lt  
A/s LT77 7044 0600 0127 1837  
AB SEB bankas  
Banko kodas 70440

### Pirkėjas

VšĮ Vilniaus universiteto ligoninė Santaros klinikos  
Santariškių g. 2, LT-08661 Vilnius  
Įmonės kodas 124364561  
PVM mok. kodas LT243645610  
A. s. LT71 7300 0100 0249 2260  
AB „Swedbank“ b. k. 73000  
Tel. (8 5) 236 5000, faks. (8 5) 236 5111

Direktorė  
Živilė Sukackienė  
A.V.




Generalinis direktorius  
Feliksas Jankevičius  
A.V.



  
Luda Orlovienė  
Viešųjų ryšių specialistė  
Tel. 8 5 270 01 01, el. paštas info@medfarmos.lt

Vyriausiasis teisininkas

  
Darius Taminskas

Viešųjų ryšių specialistas  
Tel. 8 5 270 01 01, el. paštas info@medfarmos.lt

  
Aušra Žilaitienė  
Viešųjų ryšių specialistė

Viešųjų ryšių specialistas  
Algimantas Varžgalys





20-CP-21 2020-05-14

## TECHINĖ SPECIFIKACIJA IR ĮKAINIAI

I priedas prie 2020-05- sutarties Nr.

№	Reagento arba pagalbinės priemonės pavadinimas	Specialieji reikalavimai	Gamintojo kataloginis kodas	Matavimo vienetas	Pelminiarus kieks	PVM %	1 mato vnt. įkainis Eur be PVM	1 mato vnt. įkainis Eur su PVM	Tiekėjo siūlomos pakuočės kaina, Eur be PVM	Nuoroda į gamintojo katalogą ar internetinį tinklalapį, kuriame galima patikrinti siūlomą prekęs atitikinį pirkimo dokumentuose nustatytiems reikalavimams
1	2	3	4	5	6	7	8	9	10	11
1	ZyBlack Quenching Solution	Panuoštas lipalasis autofluorescencijai mažinti.	BS-0002-8	8 ml	1	5%	45,00	47,25	47,25	<a href="https://www.zyvision.com/fr/index/zyolich/wb-0002">https://www.zyvision.com/fr/index/zyolich/wb-0002</a>
2	Fixogum Rubber Cement	Klijai tinkamys citocentiniu stikliukų panuošimui (dengiamųjų stikliukų klijavimui)	E-4005-125	125 g	21	5%	50,00	52,50	1102,5	Zr. dokumentą <a href="https://www.zyvision.com/fr/index/zyolich/wb-0002">zyvisionproductcataloguesolution/w.pdf</a> psl.: 170, 177, 232, 247
3	Fixogum Rubber Cement	Klijai tinkamys citocentiniu stikliukų panuošimui (dengiamųjų stikliukų klijavimui)	E-4005-50	50 g	1	5%	31,00	32,55	32,55	Zr. dokumentą <a href="https://www.zyvision.com/fr/index/zyolich/wb-0002">zyvisionproductcataloguesolution/w.pdf</a> psl.: 170, 177, 232, 247
4	Pepsin Solution	Panuoštas lipalasis darbu	ES-0001-1000	1000 ml	1	5%	1628,00	1709,40	1709,4	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
5	Pepsin Solution	Panuoštas lipalasis darbu/ISH	ES-0001-4	4 ml	1	5%	29,00	30,45	30,45	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
6	Pepsin Solution	Panuoštas lipalasis darbu/ISH	ES-0001-50	50 ml	1	5%	180,00	189,00	189	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
7	Pepsin Solution Set	Panuoštas lipalasis darbu/ISH	ES-0001-8	2 x 4 ml	1	5%	47,00	49,35	49,35	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
8	Cytology Pepsin Solution	Panuoštas lipalasis darbu/ISH, citologijai	ES-0002-4	4 ml	1	5%	26,00	27,30	27,3	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
9	Cytology Pepsin Solution	Panuoštas lipalasis darbu/ISH, citologijai	ES-0002-50	50 ml	1	5%	173,00	181,65	181,65	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
10	Mounting Solution (alcoholic)	Panuoštas lipalasis darbu	MT-0004-4	4 ml	1	5%	20,00	21,00	21	<a href="https://www.zyvision.com/search?searchword=MT-0004-4&amp;searchphrase=all">https://www.zyvision.com/search?searchword=MT-0004-4&amp;searchphrase=all</a>
11	DAPI/DuraTect-Solution	FISH, DAPI/Antifade	MT-0007-0.8	0,8 ml	1	5%	88,00	92,40	92,4	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
12	DAPI/DuraTect-Solution (ultra)	FISH, DAPI/Antifade	MT-0008-0.8	0,8 ml	1	5%	88,00	92,40	92,4	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
13	Heat Pretreatment Solution Citric	Heat Pretreatment Solution Citric (PT1)	PT-0001-1000	2 x 500 ml	1	5%	100,00	105,00	105	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
14	Heat Pretreatment Solution EDTA	Heat Pretreatment Solution Citric (PT2)	PT-0002-500	500 ml	1	5%	50,00	52,50	52,5	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
15	Formaldehyde Dilution Buffer Set	10x MgCl2, 50 ml; 10x PBS, 50 ml	PT-0006-100	2 x 50ml	1	5%	65,00	68,25	68,25	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
16	Wash Buffer SSC	Wash Buffer SSC	WB-0001-560	560 ml	1	5%	34,00	35,70	35,7	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
17	25x Wash Buffer A	25x Wash Buffer A	WB-0002-50	50 ml	1	5%	72,00	75,60	75,6	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
18	20x SSC Solution	20x SSC lipalasis	WB-0003-50	50 ml	1	5%	49,00	45,15	45,15	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>
19	PBS/Tween	Reagentas pipvimo etapams PBS/Tween, tab	WB-0004-1000	1 tab	1	5%	28,00	29,40	29,4	<a href="https://www.zyvision.com/search?searchword=WB-0004-1000&amp;searchphrase=all">https://www.zyvision.com/search?searchword=WB-0004-1000&amp;searchphrase=all</a>
20	20x Wash Buffer TBS	20x Wash Buffer TBS	WB-0005-50	50 ml	1	5%	47,00	49,35	49,35	<a href="https://www.zyvision.com/pro">https://www.zyvision.com/pro</a>

21	Cytology Stimpency Wash Buffer:SSC	Cytology Stimpency Wash Buffer:SSC	WB-0007-500	500 ml	1	5%	35.00	36.75	36.75	36.75	36.75	<a href="https://www.zyolight.com/br/ducts/zyolight/wb-0005">https://www.zyolight.com/br/ducts/zyolight/wb-0005</a>
22	Cytology Wash Buffer:SSC	Cytology Wash Buffer:SSC	WB-0008-500	500 ml	1	5%	35.00	36.75	36.75	36.75	36.75	<a href="https://www.zyolight.com/br/ducts/zyolight/wb-0008">https://www.zyolight.com/br/ducts/zyolight/wb-0008</a>
23	5x FlexISH Wash Buffer	5x FlexISH Wash Buffer	WB-0010-500	500 ml	1	5%	81.00	85.05	85.05	85.05	85.05	<a href="https://www.zyolight.com/br/ducts/flexish/wb-0010">https://www.zyolight.com/br/ducts/flexish/wb-0010</a>
24	Zyolight CEN 3 Probe	D3Z1 3p11-q11.1	Z-2001-200	0.2 ml	1	5%	310.00	325.50	325.50	325.50	325.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2001">https://www.zyolight.com/br/ducts/zyolight/z-2001</a>
25	Zyolight CEN 6 Probe	D6Z1 6p11.1-q11.1	Z-2002-200	0.2 ml	1	5%	310.00	325.50	325.50	325.50	325.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2002">https://www.zyolight.com/br/ducts/zyolight/z-2002</a>
26	Zyolight CEN 7 Probe	D7Z1 7q11.1	Z-2003-200	0.2 ml	1	5%	310.00	325.50	325.50	325.50	325.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2003">https://www.zyolight.com/br/ducts/zyolight/z-2003</a>
27	Zyolight CEN 8 Probe	D8Z2 8p11.1-q11.1	Z-2004-200	0.2 ml	1	5%	310.00	325.50	325.50	325.50	325.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2004">https://www.zyolight.com/br/ducts/zyolight/z-2004</a>
28	Zyolight CEN 8 Probe	D8Z2 8p11.1-q11.1	Z-2004-50	0.05 ml	1	5%	90.00	94.50	94.50	94.50	94.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2004">https://www.zyolight.com/br/ducts/zyolight/z-2004</a>
29	Zyolight CEN 11 Probe	D11Z1 11p11.1-q11.1	Z-2005-200	0.2 ml	1	5%	310.00	325.50	325.50	325.50	325.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2005">https://www.zyolight.com/br/ducts/zyolight/z-2005</a>
30	Zyolight CEN 17 Probe	D17Z1 17p11.1-q11.1	Z-2006-200	0.2 ml	1	5%	310.00	325.50	325.50	325.50	325.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2006">https://www.zyolight.com/br/ducts/zyolight/z-2006</a>
31	Zyolight CEN 18 Probe	D18Z1 18p11.1-q11.1	Z-2007-200	0.2 ml	1	5%	310.00	325.50	325.50	325.50	325.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2007">https://www.zyolight.com/br/ducts/zyolight/z-2007</a>
32	Zyolight CEN X Probe	DXZ1 Xp11.1-q11.1	Z-2008-200	0.2 ml	1	5%	310.00	325.50	325.50	325.50	325.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2008">https://www.zyolight.com/br/ducts/zyolight/z-2008</a>
33	Zyolight CEN Yq12 Probe	III DYZ1 Yq12	Z-2010-200	0.2 ml	1	5%	310.00	325.50	325.50	325.50	325.50	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2010">https://www.zyolight.com/br/ducts/zyolight/z-2010</a>
34	Zyolight SPEC MDM2/CEN 12 Dual Color Probe	MDM2 12q15 GF /CEN 12 OF: DC	Z-2013-200	0.2 ml	1	5%	960.00	1008.00	1008.00	1008.00	1008.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2013">https://www.zyolight.com/br/ducts/zyolight/z-2013</a>
35	Zyolight SPEC MDM2/CEN 12 Dual Color Probe	MDM2 12q15 GF /CEN 12 OF: DC	Z-2013-50	0.05 ml	1	5%	286.00	302.40	302.40	302.40	302.40	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2013">https://www.zyolight.com/br/ducts/zyolight/z-2013</a>
36	Zyolight SPEC MAML2 Dual Color Break Apart Probe	11q21 MAML2 GF/CEN12 (D12Z3) OF: DC BA	Z-2014-200	0.2 ml	1	5%	880.00	924.00	924.00	924.00	924.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2014">https://www.zyolight.com/br/ducts/zyolight/z-2014</a>
37	Zyolight SPEC MAML2 Dual Color Break Apart Probe	11q21 MAML2 GF/CEN12 (D12Z3) OF: DC BA	Z-2014-50	0.05 ml	1	5%	265.00	278.25	278.25	278.25	278.25	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2014">https://www.zyolight.com/br/ducts/zyolight/z-2014</a>
38	Zyolight SPEC ERBB2/CEN 17 Dual Color Probe	17q12-q21.1 ERBB2 GF/CEN 17 (D17Z1) OF: DC;	Z-2015-200	0.2 ml	1	5%	2000.00	2100.00	2100.00	2100.00	2100.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2015">https://www.zyolight.com/br/ducts/zyolight/z-2015</a>
39	Zyolight SPEC ERBB2/CEN 17 Dual Color Probe	17q12-q21.1 ERBB2 GF/CEN 17 (D17Z1) OF: DC;	Z-2015-50	0.05 ml	1	5%	560.00	588.00	588.00	588.00	588.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2015">https://www.zyolight.com/br/ducts/zyolight/z-2015</a>
40	Zyolight CEN X7q12 Dual Color Probe	DXZ1/III DYZ1 Xp11.1-q11.1/7q12	Z-2016-200	0.2 ml	1	5%	300.00	315.00	315.00	315.00	315.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2016">https://www.zyolight.com/br/ducts/zyolight/z-2016</a>
41	Zyolight CEN X7q12 Dual Color Probe	DXZ1/III DYZ1 Xp11.1-q11.1/7q12	Z-2016-50	0.05 ml	1	5%	100.00	105.00	105.00	105.00	105.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2016">https://www.zyolight.com/br/ducts/zyolight/z-2016</a>
42	Zyolight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe	13q14.11 FOXO1 OF/p36.1/PAX3 GF: DC; RMS I Probe SPEC (t2:13): DC FP	Z-2018-200	0.2 ml	1	5%	880.00	924.00	924.00	924.00	924.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2018">https://www.zyolight.com/br/ducts/zyolight/z-2018</a>
43	Zyolight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe	13q14.11 FOXO1 OF/p36.1/PAX3 GF: DC; RMS I Probe SPEC (t2:13): DC FP	Z-2018-50	0.05 ml	1	5%	264.00	277.20	277.20	277.20	277.20	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2018">https://www.zyolight.com/br/ducts/zyolight/z-2018</a>
44	Zyolight SPEC FOXO1/PAX7 Dual Color Single Fusion Probe	13q14.11 FOXO1 OF/p36.13 PAX7 GF: DC; RMS II Probe SPEC (t1:13): DC FP	Z-2019-200	0.2 ml	1	5%	880.00	924.00	924.00	924.00	924.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2019">https://www.zyolight.com/br/ducts/zyolight/z-2019</a>
45	Zyolight SPEC FOXO1/PAX7 Dual Color Single Fusion Probe	13q14.11 FOXO1 OF/p36.13 PAX7 GF: DC; RMS II Probe SPEC (t1:13): DC FP	Z-2019-50	0.05 ml	1	5%	264.00	277.20	277.20	277.20	277.20	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2019">https://www.zyolight.com/br/ducts/zyolight/z-2019</a>
46	Zyolight SPEC ERBB2/CEN 17 Dual Color Probe Kit	17q12-q21.1 ERBB2 GF/CEN 17 (D17Z1) OF: DC kitkiny	Z-2020-20	20 testq	1	5%	2440.00	2562.00	2562.00	2562.00	2562.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2020">https://www.zyolight.com/br/ducts/zyolight/z-2020</a>
47	Zyolight SPEC ERBB2/CEN 17 Dual Color Probe Kit	17q12-q21.1 ERBB2 GF/CEN 17 (D17Z1) OF: DC kitkiny	Z-2020-5	5 testai	1	5%	660.00	693.00	693.00	693.00	693.00	<a href="https://www.zyolight.com/br/ducts/zyolight/z-2020">https://www.zyolight.com/br/ducts/zyolight/z-2020</a>

48	ZyeloLight FISH-Tissue Implementation Kit	Heat Pretreatment Solution Clinic, Pepsin Solution, Wash Buffer SSC, 25x Wash Buffer A, DAPI/DuraTect-Solution	Z-2028-20	20 test/s	1	5%	324.00	340.20	340.2	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2028">https://www.zylovision.com/pro-ducts/zyolight/z-2028</a>
49	ZyeloLight FISH-Tissue Implementation Kit	Heat Pretreatment Solution Clinic, Pepsin Solution, Wash Buffer SSC, 25x Wash Buffer A, DAPI/DuraTect-Solution	Z-2028-5	5 test/s	1	5%	100.00	105.00	105	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2028">https://www.zylovision.com/pro-ducts/zyolight/z-2028</a>
50	ZyeloLight SPEC EGFR/CEN 7 Dual Color Probe	7p11.2 EGFRCEN 7 (D7Z1) OF, DC	Z-2033-200	0.2 ml	1	5%	985.00	1034.25	1034.25	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2033">https://www.zylovision.com/pro-ducts/zyolight/z-2033</a>
51	ZyeloLight SPEC EGFR/CEN 7 Dual Color Probe	7p11.2 EGFRCEN 7 (D7Z1) OF, DC	Z-2033-50	0.05 ml	1	5%	305.00	320.25	320.25	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2033">https://www.zylovision.com/pro-ducts/zyolight/z-2033</a>
52	ZyeloLight SPEC 2q11 Probe	2q11.2	Z-2049-200	0.2 ml	1	5%	305.00	320.25	320.25	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2049">https://www.zylovision.com/pro-ducts/zyolight/z-2049</a>
53	ZyeloLight CEN 12 Probe	D12Z3 12p11.1-q11	Z-2050-200	0.2 ml	1	5%	250.00	262.50	262.5	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2050">https://www.zylovision.com/pro-ducts/zyolight/z-2050</a>
54	ZyeloLight SPEC ERBB3/CEN 12 Dual Color Probe	12q13.2-q13.3 ERBB3 OF CEN12 (D12Z3) GF, DC	Z-2056-200	0.2 ml	1	5%	1100.00	1155.00	1155	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2056">https://www.zylovision.com/pro-ducts/zyolight/z-2056</a>
55	ZyeloLight SPEC ERBB4/2q11 Dual Color Probe	2q34 ERBB4 GF/2q11 OF, AFF3 (AF4/FMR2) 2q11.2, DC	Z-2057-200	0.2 ml	1	5%	1100.00	1155.00	1155	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2057">https://www.zylovision.com/pro-ducts/zyolight/z-2057</a>
56	ZyeloLight SPEC FHT/CEN 3 Dual Color Probe	3p14.2 FHT GF/CEN 3 (D3Z1) OF, DC	Z-2062-200	0.2 ml	1	5%	1100.00	1155.00	1155	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2062">https://www.zylovision.com/pro-ducts/zyolight/z-2062</a>
57	ZyeloLight SPEC CDKN2A/CEN 9 Dual Color Probe	9p21.3 CDKN2A GF/CEN 9 (D9Z3) OF, DC	Z-2063-200	0.2 ml	1	5%	1000.00	1050.00	1050	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2063">https://www.zylovision.com/pro-ducts/zyolight/z-2063</a>
58	ZyeloLight SPEC CDKN2A/CEN 9 Dual Color Probe	9p21.3 CDKN2A GF/CEN 9 (D9Z3) OF, DC	Z-2063-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2063">https://www.zylovision.com/pro-ducts/zyolight/z-2063</a>
59	ZyeloLight CEN 9 Probe	11l D9Z3 9q12	Z-2067-200	0.2 ml	1	5%	250.00	262.50	262.5	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2067">https://www.zylovision.com/pro-ducts/zyolight/z-2067</a>
60	ZyeloLight SPEC ESR1/CEN 6 Dual Color Probe	6q25.1 ESR1 GF/CEN 6 (D6Z1) OF, DC	Z-2069-200	0.2 ml	1	5%	900.00	945.00	945	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2069">https://www.zylovision.com/pro-ducts/zyolight/z-2069</a>
61	ZyeloLight SPEC ESR1/CEN 6 Dual Color Probe	6q25.1 ESR1 GF/CEN 6 (D6Z1) OF, DC	Z-2069-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2069">https://www.zylovision.com/pro-ducts/zyolight/z-2069</a>
62	ZyeloLight SPEC CCND1/CEN 11 Dual Color Probe	11q13.3 CCND1 GF/CEN 11 (D11Z1) OF, DC	Z-2071-200	0.2 ml	1	5%	750.00	787.50	787.5	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2071">https://www.zylovision.com/pro-ducts/zyolight/z-2071</a>
63	ZyeloLight SPEC CCND1/CEN 11 Dual Color Probe	11q13.3 CCND1 GF/CEN 11 (D11Z1) OF, DC	Z-2071-50	0.05 ml	1	5%	220.00	231.00	231	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2071">https://www.zylovision.com/pro-ducts/zyolight/z-2071</a>
64	ZyeloLight SPEC FGFR1/CEN 8 Dual Color Probe	8p11.23-p11.22FGFR1 GF/CEN 8(D8Z2) OF, DC	Z-2072-200	0.2 ml	1	5%	900.00	945.00	945	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2072">https://www.zylovision.com/pro-ducts/zyolight/z-2072</a>
65	ZyeloLight SPEC FGFR1/CEN 8 Dual Color Probe	8p11.23-p11.22FGFR1 GF/CEN 8(D8Z2) OF, DC	Z-2072-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2072">https://www.zylovision.com/pro-ducts/zyolight/z-2072</a>
66	ZyeloLight SPEC MYCN/2q11 Dual Color Probe	2q24.3 MYCN GF/2q11 AFF3 (AF4/FMR2) 2q11.2 OF, DC	Z-2074-200	0.2 ml	2	5%	1200.000	1260.00	2520	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2074">https://www.zylovision.com/pro-ducts/zyolight/z-2074</a>
67	ZyeloLight SPEC MYCN/2q11 Dual Color Probe	2q24.3 MYCN GF/2q11 AFF3 (AF4/FMR2) 2q11.2 OF, DC	Z-2074-50	0.05 ml	1	5%	340.000	357.00	357	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2074">https://www.zylovision.com/pro-ducts/zyolight/z-2074</a>
68	ZyeloLight SPEC 1p36/1q25 Dual Color Probe	1p36 1p36.31 OF/1q25 1q25.3 GF, DC	Z-2075-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2075">https://www.zylovision.com/pro-ducts/zyolight/z-2075</a>
69	ZyeloLight SPEC 1p36/1q25 Dual Color Probe	1p36 1p36.31 OF/1q25 1q25.3 GF, DC	Z-2075-50	0.05 ml	1	5%	275.00	288.75	288.75	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2075">https://www.zylovision.com/pro-ducts/zyolight/z-2075</a>
70	ZyeloLight SPEC 19q13/19p13 Dual Color Probe	19q13 19q13.32-q13.33 OF/19p13 19p13.3 GF, DC	Z-2076-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2076">https://www.zylovision.com/pro-ducts/zyolight/z-2076</a>
71	ZyeloLight SPEC 19q13/19p13 Dual Color Probe	19q13 19q13.32-q13.33 OF/19p13 19p13.3 GF, DC	Z-2076-50	0.05 ml	1	5%	275.00	288.75	288.75	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2076">https://www.zylovision.com/pro-ducts/zyolight/z-2076</a>
72	ZyeloLight CEN 17/SPEC ERBB2 Dual Color Probe	CEN 17/SPEC ERBB2 DC	Z-2077-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2077">https://www.zylovision.com/pro-ducts/zyolight/z-2077</a>
73	ZyeloLight CEN 17/SPEC ERBB2 Dual Color Probe	CEN 17/SPEC ERBB2 DC	Z-2077-50	0.05 ml	1	5%	275.00	288.75	288.75	<a href="https://www.zylovision.com/pro-ducts/zyolight/z-2077">https://www.zylovision.com/pro-ducts/zyolight/z-2077</a>
74	ZyeloLight SPEC PTEN/CEN 10 Dual Color Probe	10q23.2-q23.31 PTEN GF/CEN 10 (D10Z1) OF, DC	Z-2078-200	0.2 ml	1	5%	780.00	819.00	819	<a href="https://www.zylovision.com/pro-ducts/zyolight/sz-2078">https://www.zylovision.com/pro-ducts/zyolight/sz-2078</a>

75	ZyLight SPEC PTEN/CEN 10 Dual Color Probe	10q23.2-q23.31 PTEN GF/CEN 10 (D10Z1) OF, DC	Z-2078-50	0.05 ml	1	5%	240.00	252.00	252	<a href="https://www.zyovision.com/products/zyolight/z-2078">https://www.zyovision.com/products/zyolight/z-2078</a>
76	ZyLight CEN 10 Probe	D10Z1 10p11.1-q11.1	Z-2079-200	0.2 ml	1	5%	250.00	262.50	262.5	<a href="https://www.zyovision.com/products/zyolight/z-2079">https://www.zyovision.com/products/zyolight/z-2079</a>
77	ZyLight SPEC MDM4/1p12 Dual Color Probe	1q32.1MDM4 GF/1p12 OF, DC	Z-2080-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zyovision.com/products/zyolight/z-2080">https://www.zyovision.com/products/zyolight/z-2080</a>
78	ZyLight SPEC CDKN2A/CEN 3/7/17 Quaduple Color Probe	9p21.3 CDKN2A GF/CEN 3 (D3Z1) RF/ CEN 7 (D7Z1) GF/ CEN 17 (D17Z1) BF, QC	Z-2081-200	0.2 ml	1	5%	2300.00	2415.00	2415	<a href="https://www.zyovision.com/products/zyolight/z-2081">https://www.zyovision.com/products/zyolight/z-2081</a>
79	ZyLight SPEC CDKN2A/CEN 3/7/17 Quaduple Color Probe	9p21.3 CDKN2A GF/CEN 3 (D3Z1) RF/ CEN 7 (D7Z1) GF/ CEN 17 (D17Z1) BF, QC	Z-2081-50	0.05 ml	1	5%	690.00	724.50	724.5	<a href="https://www.zyovision.com/products/zyolight/z-2081">https://www.zyovision.com/products/zyolight/z-2081</a>
80	ZyLight SPEC FGFR3/4p11 Dual Color Probe	4p16.3 FGFR3 GF/4p11 ZR1 OF, DC	Z-2082-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zyovision.com/products/zyolight/z-2082">https://www.zyovision.com/products/zyolight/z-2082</a>
81	ZyLight SPEC 4p11 Probe	4p11	Z-2083-200	0.2 ml	1	5%	290.00	304.50	304.5	<a href="https://www.zyovision.com/products/zyolight/z-2083">https://www.zyovision.com/products/zyolight/z-2083</a>
82	ZyLight SPEC VHL/CEN 3 Dual Color Probe	3p25.3 VHL GF/CEN 3 (D3Z1) OF, DC	Z-2084-200	0.2 ml	1	5%	950.00	997.50	997.5	<a href="https://www.zyovision.com/products/zyolight/z-2084">https://www.zyovision.com/products/zyolight/z-2084</a>
83	ZyLight SPEC 13q12 Probe	13q12.11	Z-2085-200	0.2 ml	1	5%	290.00	304.50	304.5	<a href="https://www.zyovision.com/products/zyolight/z-2085">https://www.zyovision.com/products/zyolight/z-2085</a>
84	ZyLight SPEC 21q22 Probe	21q22.13-q22.2	Z-2086-200	0.2 ml	1	5%	290.00	304.50	304.5	<a href="https://www.zyovision.com/products/zyolight/z-2086">https://www.zyovision.com/products/zyolight/z-2086</a>
85	ZyLight SPEC MET/CEN 7 Dual Color Probe	7q31.2 MET GF/CEN (D7Z1) 7 OF, DC	Z-2087-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zyovision.com/products/zyolight/z-2087">https://www.zyovision.com/products/zyolight/z-2087</a>
86	ZyLight SPEC MET/CEN 7 Dual Color Probe	7q31.2 MET GF/CEN (D7Z1) 7 OF, DC	Z-2087-50	0.05 ml	1	5%	290.00	304.50	304.5	<a href="https://www.zyovision.com/products/zyolight/z-2087">https://www.zyovision.com/products/zyolight/z-2087</a>
87	ZyLight SPEC MYC Dual Color Break Apart Probe	MYC DC BA: 8q24.21 proximal OF and distal GF	Z-2089-200	0.2 ml	5	5%	1300.000	1365.00	6825	<a href="https://www.zyovision.com/products/zyolight/z-2089">https://www.zyovision.com/products/zyolight/z-2089</a>
88	ZyLight SPEC MYC Dual Color Break Apart Probe	MYC DC BA: 8q24.21 proximal OF and distal GF	Z-2090-50	0.05 ml	1	5%	345.000	362.25	362.25	<a href="https://www.zyovision.com/products/zyolight/z-2090">https://www.zyovision.com/products/zyolight/z-2090</a>
89	ZyLight SPEC TERT/5q31 Dual Color Probe	5p15.33 TERT GF/5q31 OF, DC	Z-2091-200	0.2 ml	1	5%	850.00	892.50	892.5	<a href="https://www.zyovision.com/products/zyolight/z-2091">https://www.zyovision.com/products/zyolight/z-2091</a>
90	ZyLight SPEC TERT/5q31 Dual Color Probe	5p15.33 TERT GF/5q31 OF, DC	Z-2091-50	0.05 ml	1	5%	250.00	262.50	262.5	<a href="https://www.zyovision.com/products/zyolight/z-2091">https://www.zyovision.com/products/zyolight/z-2091</a>
91	ZyLight SPEC MYC/CEN 8 Dual Color Probe	8q24.21 MYC GF/CEN 8 (D8Z2) OF, DC	Z-2092-200	0.2 ml	1	5%	750.00	787.50	787.5	<a href="https://www.zyovision.com/products/zyolight/z-2092">https://www.zyovision.com/products/zyolight/z-2092</a>
92	ZyLight SPEC MYC/CEN 8 Dual Color Probe	8q24.21 MYC GF/CEN 8 (D8Z2) OF, DC	Z-2092-50	0.05 ml	1	5%	220.00	231.00	231	<a href="https://www.zyovision.com/products/zyolight/z-2092">https://www.zyovision.com/products/zyolight/z-2092</a>
93	ZyLight SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe	17q12-q21.1 ERBB2 GF/17q21.2 TOP2A OF/CEN 17 (D17Z1) BF, TIC	Z-2093-200	0.2 ml	1	5%	1500.00	1575.00	1575	<a href="https://www.zyovision.com/products/zyolight/z-2093">https://www.zyovision.com/products/zyolight/z-2093</a>
94	ZyLight SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe	17q12-q21.1 ERBB2 GF/17q21.2 TOP2A OF/CEN 17 (D17Z1) BF, TIC	Z-2093-50	0.05 ml	1	5%	500.00	525.00	525	<a href="https://www.zyovision.com/products/zyolight/z-2093">https://www.zyovision.com/products/zyolight/z-2093</a>
95	ZyLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe	13q12.11/18p11.1/21q22.13-q22.2 Green/Blue/Orange	Z-2095-200	0.2 ml	1	5%	900.00	945.00	945	<a href="https://www.zyovision.com/products/zyolight/z-2095">https://www.zyovision.com/products/zyolight/z-2095</a>
96	ZyLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe	13q12.11/18p11.1/21q22.13-q22.2 Green/Blue/Orange	Z-2095-50	0.05 ml	1	5%	250.00	262.50	262.5	<a href="https://www.zyovision.com/products/zyolight/z-2095">https://www.zyovision.com/products/zyolight/z-2095</a>
97	ZyLight SPEC EWSR1 Dual Color Break Apart Probe	EWSR1 DC 22q12.1-q12.2 EWSR1 OF proximal and extends inward into intron 4 /GF distal	Z-2096-200	0.2 ml	1	5%	1015.00	1065.75	1065.75	<a href="https://www.zyovision.com/products/zyolight/z-2096">https://www.zyovision.com/products/zyolight/z-2096</a>
98	ZyLight SPEC EWSR1 Dual Color Break Apart Probe	EWSR1 DC 22q12.1-q12.2 EWSR1 OF proximal and extends inward into intron 4 /GF distal	Z-2096-50	0.05 ml	1	5%	330.00	346.50	346.5	<a href="https://www.zyovision.com/products/zyolight/z-2096">https://www.zyovision.com/products/zyolight/z-2096</a>

99	ZyoLight SPEC SS18 Dual Color Break Apart Probe	SS18 DC BA 18q11.2: OF distal to the SS18 gene/GF proximal to that gene	Z-2097-200	0.2 ml	1	5%	820.00	861.00	861	<a href="https://www.zyovision.com/products/zyolight/z-2097">https://www.zyovision.com/products/zyolight/z-2097</a>
100	ZyoLight SPEC SS18 Dual Color Break Apart Probe	SS18 DC BA 18q11.2: OF distal to the SS18 gene/GF proximal to that gene	Z-2097-50	0.05 ml	1	5%	250.00	262.50	262,5	<a href="https://www.zyovision.com/products/zyolight/z-2097">https://www.zyovision.com/products/zyolight/z-2097</a>
101	ZyoLight FISH-Cytology Implementation Kit	Cytology Pepsin Solution, 20x Wash Buffer TBS, 10x MgCl <sub>2</sub> , 10x PBS, Cytology Stingency Wash Buffer SSC, Cytology Wash Buffer SSC, DAPI/DuraTect-Solution	Z-2099-20	20 lyriml	1	5%	324.00	340.20	340,2	<a href="https://www.zyovision.com/products/zyolight/z-2099">https://www.zyovision.com/products/zyolight/z-2099</a>
102	ZyoLight SPEC DDIT3 Dual Color Break Apart Probe	12q13.3-q14.1 DDIT3 DC BA: OF proximal to the DDIT3 gene/GF distal to that gene.	Z-2100-200	0.2 ml	1	5%	725.00	761.25	761,25	<a href="https://www.zyovision.com/products/zyolight/z-2100">https://www.zyovision.com/products/zyolight/z-2100</a>
103	ZyoLight SPEC DDIT3 Dual Color Break Apart Probe	12q13.3-q14.1 DDIT3 DC BA: OF proximal to the DDIT3 gene/GF distal to that gene.	Z-2100-50	0.05 ml	1	5%	220.00	231.00	231	<a href="https://www.zyovision.com/products/zyolight/z-2100">https://www.zyovision.com/products/zyolight/z-2100</a>
104	ZyoLight SPEC 1p12 Probe	1p12	Z-2101-200	0.2 ml	1	5%	295.00	309.75	309,75	<a href="https://www.zyovision.com/products/zyolight/z-2101">https://www.zyovision.com/products/zyolight/z-2101</a>
105	ZyoLight SPEC VHL/1p12/CEN 7/12 Quadruple Color Probe	3p25.3 VHL GF/1p12 RF/CEN 7 (D7Z1), G6F/CEN 17 (D17Z1) BF, QC	Z-2102-200	0.2 ml	1	5%	2100.00	2205.00	2205	<a href="https://www.zyovision.com/products/zyolight/z-2102">https://www.zyovision.com/products/zyolight/z-2102</a>
106	ZyoLight SPEC CDK4/CEN 12 Dual Color Probe	12q13.3-q14.1 CDK4 GF/CEN 12 (D12Z3) OF, DC	Z-2103-200	0.2 ml	1	5%	900.00	945.00	945	<a href="https://www.zyovision.com/products/zyolight/z-2103">https://www.zyovision.com/products/zyolight/z-2103</a>
107	ZyoLight SPEC CDK4/CEN 12 Dual Color Probe	12q13.3-q14.1 CDK4 GF/CEN 12 (D12Z3) OF, DC	Z-2103-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zyovision.com/products/zyolight/z-2103">https://www.zyovision.com/products/zyolight/z-2103</a>
108	ZyoLight Aneuploidy Panel XY and 13/18/21	Incl. ZyoLight CEN X/Y/q12 Dual Color Probe, 0.2 ml (Z-2016-200); ZyoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe, 0.2 ml (Z-2095-200) Panel XY and 13/18/21	Z-2104-20	2x 0.2 ml	1	5%	1200.00	1260.00	1260	<a href="https://www.zyovision.com/products/zyolight/gene-index">https://www.zyovision.com/products/zyolight/gene-index</a>
109	ZyoLight Aneuploidy Panel XY and 13/18/21	Incl. ZyoLight CEN X/Y/q12 Dual Color Probe, 0.2 ml (Z-2016-200); ZyoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe, 0.2 ml (Z-2095-200) Panel XY and 13/18/21	Z-2104-5	2x 0.05 ml	1	5%	320.00	336.00	336	<a href="https://www.zyovision.com/products/zyolight/gene-index">https://www.zyovision.com/products/zyolight/gene-index</a>
110	ZyoLight SPEC MYC/IGH Dual Color Dual Fusion Probe	MYC OF/IGH GF, DC	Z-2105-200	0.2 ml	1	5%	955.00	1002.75	1002,75	<a href="https://www.zyovision.com/products/zyolight/z-2105">https://www.zyovision.com/products/zyolight/z-2105</a>
111	ZyoLight SPEC MYC/IGH Dual Color Dual Fusion Probe	MYC OF/IGH GF, DC	Z-2105-50	0.05 ml	1	5%	300.00	315.00	315	<a href="https://www.zyovision.com/products/zyolight/z-2105">https://www.zyovision.com/products/zyolight/z-2105</a>
112	ZyoLight SPEC EGR1/5p15 Dual Color Probe	5q31.2 EGR1 GF/5p15 5p15.33 OF, DC	Z-2107-200	0.2 ml	1	5%	730.00	766.50	766,5	<a href="https://www.zyovision.com/products/zyolight/z-2107">https://www.zyovision.com/products/zyolight/z-2107</a>
113	ZyoLight SPEC EGR1/5p15 Dual Color Probe	5q31.2 EGR1 GF/5p15 5p15.33 OF, DC	Z-2107-50	0.05 ml	1	5%	230.00	241.50	241,5	<a href="https://www.zyovision.com/products/zyolight/z-2107">https://www.zyovision.com/products/zyolight/z-2107</a>
114	ZyoLight SPEC CCND1 Dual Color Break Apart Probe	CCND1 11q13.2-q13.3 DC BA: OF distal to and covers the CCND1 gene/GF proximal to that gene.	Z-2108-200	0.2 ml	2	5%	1065.000	1118.25	3354,75	<a href="https://www.zyovision.com/products/zyolight/z-2108">https://www.zyovision.com/products/zyolight/z-2108</a>
115	ZyoLight SPEC CCND1 Dual Color Break Apart Probe	CCND1 11q13.2-q13.3 DC BA: OF distal to and covers the CCND1 gene/GF proximal to that gene.	Z-2108-50	0.05 ml	1	5%	330.000	346.50	346,5	<a href="https://www.zyovision.com/products/zyolight/z-2108">https://www.zyovision.com/products/zyolight/z-2108</a>
116	ZyoLight SPEC TFE3 Dual Color Break Apart Probe	Xp11.23 TFE3 DC BA: OF probe hybridizes distal to the TFE3 gene/GF proximal to that gene.	Z-2109-200	0.2 ml	2	5%	1150.000	1207.50	2415	<a href="https://www.zyovision.com/products/zyolight/z-2109">https://www.zyovision.com/products/zyolight/z-2109</a>
117	ZyoLight SPEC TFE3 Dual Color Break Apart Probe	Xp11.23 TFE3 DC BA: OF probe hybridizes distal to the TFE3 gene/GF proximal to that gene.	Z-2109-50	0.05 ml	1	5%	330.000	346.50	346,5	<a href="https://www.zyovision.com/products/zyolight/z-2109">https://www.zyovision.com/products/zyolight/z-2109</a>
118	ZyoLight SPEC IGH Dual Color Break Apart Probe	14q32.33 IGH DC BA: OF proximal/GF distal to the constant regions of the IGH locus	Z-2110-200	0.2 ml	1	5%	730.00	766.50	766,5	<a href="https://www.zyovision.com/products/zyolight/z-2110">https://www.zyovision.com/products/zyolight/z-2110</a>
119	ZyoLight SPEC IGH Dual Color Break Apart Probe	14q32.33 IGH DC BA: OF proximal/GF distal to the constant regions of the IGH locus	Z-2110-50	0.05 ml	1	5%	230.00	241.50	241,5	<a href="https://www.zyovision.com/products/zyolight/z-2110">https://www.zyovision.com/products/zyolight/z-2110</a>
120	ZyoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe	BCR GF/ABL1 OF, DC	Z-2111-200	0.2 ml	1	5%	800.00	840.00	840	<a href="https://www.zyovision.com/products/zyolight/z-2111">https://www.zyovision.com/products/zyolight/z-2111</a>
121	ZyoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe	BCR GF/ABL1 OF, DC	Z-2111-50	0.05 ml	1	5%	240.00	252.00	252	<a href="https://www.zyovision.com/products/zyolight/z-2111">https://www.zyovision.com/products/zyolight/z-2111</a>

122	ZyloLight SPEC RUNX1/RUNX1T1 Dual Color Fusion Probe	RUNX1 GF/RUNX1T1 OF, DC DF.	Z-2112-200	0.2 ml	1	5%	800,00	840,00	840	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2112">https://www.zylovision.com/probe/ducts/zytolight/z-2112</a>
123	ZyloLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe	RUNX1 GF/RUNX1T1 OF, DC DF.	Z-2112-50	0.05 ml	1	5%	250,00	262,50	262,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2112">https://www.zylovision.com/probe/ducts/zytolight/z-2112</a>
124	ZyloLight SPEC PML/RARA Dual Color Dual Fusion Probe	PML OF/RARA GF, DC DF.	Z-2113-200	0.2 ml	1	5%	770,00	808,50	808,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2113">https://www.zylovision.com/probe/ducts/zytolight/z-2113</a>
125	ZyloLight SPEC PML/RARA Dual Color Dual Fusion Probe	PML OF/RARA GF, DC DF.	Z-2113-50	0.05 ml	1	5%	250,00	262,50	262,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2113">https://www.zylovision.com/probe/ducts/zytolight/z-2113</a>
126	ZyloLight SPEC BCL2/IGH Dual Color Dual Fusion Probe	BCL2 OF/IGH GF, DC DF	Z-2114-200	0.2 ml	1	5%	990,00	1039,50	1039,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2114">https://www.zylovision.com/probe/ducts/zytolight/z-2114</a>
127	ZyloLight SPEC BCL2/IGH Dual Color Dual Fusion Probe	BCL2 OF/IGH GF, DC DF	Z-2114-50	0.05 ml	1	5%	350,00	367,50	367,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2114">https://www.zylovision.com/probe/ducts/zytolight/z-2114</a>
128	ZyloLight SPEC KRAS/CEN 12 Dual Color Probe	12p12.1 KRAS GF/CEN 12 (D1223) OF, DC	Z-2115-200	0.2 ml	1	5%	930,00	976,50	976,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2115">https://www.zylovision.com/probe/ducts/zytolight/z-2115</a>
129	ZyloLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe	BCR GF/ABL1 OF, DC	Z-2116-200	0.2 ml	1	5%	875,00	918,75	918,75	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2116">https://www.zylovision.com/probe/ducts/zytolight/z-2116</a>
130	ZyloLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe	BCR GF/ABL1 OF, DC	Z-2116-50	0.05 ml	1	5%	263,00	276,15	276,15	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2116">https://www.zylovision.com/probe/ducts/zytolight/z-2116</a>
131	ZyloLight SPEC ALK/EML4 TriCheck Probe	2p23.2 ALK OF GF 2p23.1-p23.2 EML4 2p21 BF, TrC	Z-2117-200	0.2 ml	1	5%	1300,00	1365,00	1365	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2117">https://www.zylovision.com/probe/ducts/zytolight/z-2117</a>
132	ZyloLight SPEC ALK/EML4 TriCheck Probe	2p23.2 ALK OF GF 2p23.1-p23.2 EML4 2p21 BF, TrC	Z-2117-50	0.05 ml	1	5%	384,00	403,20	403,2	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2117">https://www.zylovision.com/probe/ducts/zytolight/z-2117</a>
133	ZyloLight SPEC CCND1 Break Apart/2q11/CEN 6 Quadruple Color Probe	COND1 Break Apart/2q11/CEN 6 Quadruple Color	Z-2118-200	0.2 ml	1	5%	2240,00	2352,00	2352	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2118">https://www.zylovision.com/probe/ducts/zytolight/z-2118</a>
134	ZyloLight SPEC PDGFB Dual Color Break Apart Probe	PDGFB 22q13.1 DC BA: OF/GF	Z-2119-200	0.2 ml	1	5%	835,00	876,75	876,75	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2119">https://www.zylovision.com/probe/ducts/zytolight/z-2119</a>
135	ZyloLight SPEC PDGFB Dual Color Break Apart Probe	PDGFB 22q13.1 DC BA: OF/GF	Z-2119-50	0.05 ml	1	5%	250,00	262,50	262,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2119">https://www.zylovision.com/probe/ducts/zytolight/z-2119</a>
136	ZyloLight CEN XY Dual Color Probe	DXZ1/DY23 Xp11.1-q11.1/Yp11.1-q11.1	Z-2120-200	0.2 ml	1	5%	350,00	367,50	367,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2120">https://www.zylovision.com/probe/ducts/zytolight/z-2120</a>
137	ZyloLight SPEC COL1A1 Dual Color Break Apart Probe	17q21.33 COL1A1 DC BA: OF/GF	Z-2121-200	0.2 ml	1	5%	820,00	861,00	861	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2121">https://www.zylovision.com/probe/ducts/zytolight/z-2121</a>
138	ZyloLight SPEC FGF2/CEN 10 Dual Color Probe	10q26.12-q26.13 FGF2 GF/CEN 10 (D1021) OF, DC.	Z-2122-200	0.2 ml	1	5%	950,00	997,50	997,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2122">https://www.zylovision.com/probe/ducts/zytolight/z-2122</a>
139	ZyloLight CEN Y (DY23) Probe	DY23 Yp11.1-q11.1	Z-2123-200	0.2 ml	1	5%	250,00	262,50	262,5	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2123">https://www.zylovision.com/probe/ducts/zytolight/z-2123</a>
140	ZyloLight SPEC ALK Dual Color Break Apart Probe	2p23.1-p23.2 ALK DC BA: OF/GF	Z-2124-200	0.2 ml	1	5%	3020,00	3171,00	3171	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2124">https://www.zylovision.com/probe/ducts/zytolight/z-2124</a>
141	ZyloLight SPEC ALK Dual Color Break Apart Probe	2p23.1-p23.2 ALK DC BA: OF/GF	Z-2124-50	0.05 ml	1	5%	800,00	840,00	840	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2124">https://www.zylovision.com/probe/ducts/zytolight/z-2124</a>
142	ZyloLight SPEC CCND1/IGH Dual Color Dual Fusion Probe	CCND1 OF/IGH GF, DC DF	Z-2125-200	0.2 ml	1	5%	722,00	758,10	758,1	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2125">https://www.zylovision.com/probe/ducts/zytolight/z-2125</a>
143	ZyloLight SPEC CCND1/IGH Dual Color Dual Fusion Probe	CCND1 OF/IGH GF, DC DF	Z-2125-50	0.05 ml	1	5%	220,00	231,00	231	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2125">https://www.zylovision.com/probe/ducts/zytolight/z-2125</a>
144	ZyloLight SPEC SOX2/CEN 3 Dual Color Probe	3q26.33 SOX2 GF/CEN 3 (D321) OF, DC	Z-2127-200	0.2 ml	1	5%	920,00	966,00	966	<a href="https://www.zylovision.com/probe/ducts/zytolight/z-2127">https://www.zylovision.com/probe/ducts/zytolight/z-2127</a>

145	ZyeloLight SPEC FUS Dual Color Break Apart Probe	16p11.2 FUS DC BA: OF/GF	Z-2150-50	0.05 ml	1	5%	220.00	231.00	231	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2130</a>
146	ZyeloLight SPEC KIF5B Dual Color Break Apart Probe	10p11.22 KIF5B DC BA: OF/GF	Z-2151-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2131</a>
147	ZyeloLight SPEC JAZF1 Dual Color Break Apart Probe	7p15.1-p15.2 JAZF1 DC BA: OF/GF	Z-2152-50	0.05 ml	8	5%	367.00	385.35	3082.8	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2132</a>
148	ZyeloLight SPEC ERG/MPRSS2 Trichetec: OF/GF/BF Probe	21q22.13-q22.3 ERG/MPRSS2 Trichetec: OF/GF/BF	Z-2155-200	0.2 ml	1	5%	1280.00	1344.00	1344	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2135</a>
149	ZyeloLight SPEC EML4 Dual Color Break Apart Probe	2p21 EML4 DC BA: OF/GF	Z-2156-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2136</a>
150	ZyeloLight SPEC CARS Dual Color Break Apart Probe	11p15.4-p15.5 CARS DC BA: OF/GF	Z-2157-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2137</a>
151	ZyeloLight SPEC ERG Dual Color Break Apart Probe	21q22.13-q22.2 ERG DC BA: OF/GF	Z-2158-200	0.2 ml	1	5%	930.00	976.50	976.5	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2138</a>
152	ZyeloLight SPEC FOXO1 Dual Color Break Apart Probe	13q14.11 FOXO1 DC BA: OF/GF	Z-2159-50	0.05 ml	1	5%	255.00	267.75	267.75	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2139</a>
153	ZyeloLight SPEC PIK3CA/CEN 3 Dual Color Probe	3q26.32-q26.33 PIK3CA GF /CEN 3 (D321) OF: DC	Z-2140-200	0.2 ml	1	5%	925.00	971.25	971.25	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2140</a>
154	ZyeloLight SPEC WT1 Dual Color Break Apart Probe	11p13 WT1 DC BA: OF/GF	Z-2142-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2142</a>
155	ZyeloLight SPEC MYB Dual Color Break Apart Probe	6q23.2-q23.3 MYB DC BA: OF/GF	Z-2143-200	0.2 ml	1	5%	880.00	924.00	924	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2143</a>
156	ZyeloLight SPEC MYB Dual Color Break Apart Probe	6q23.2-q23.3 MYB DC BA: OF/GF	Z-2143-50	0.05 ml	1	5%	270.00	283.50	283.5	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2143</a>
157	ZyeloLight SPEC ROS1 Dual Color Break Apart Probe	6q22.1 ROS1 DC BA: OF/GF	Z-2144-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2144</a>
158	ZyeloLight SPEC ROS1 Dual Color Break Apart Probe	6q22.1 ROS1 DC BA: OF/GF	Z-2144-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2144</a>
159	ZyeloLight SPEC NR4A3 Dual Color Break Apart Probe	9q22.33-q31.1 NR4A3 DC BA: OF/GF	Z-2145-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2145</a>
160	ZyeloLight SPEC BIRC3/MAL T1 Dual Color Fusion Probe	11q22.1-q22.2 BIRC3 GF /1q21.31-q21.32 MAL T1 OF: DCD	Z-2146-200	0.2 ml	1	5%	810.00	850.50	850.5	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2146</a>
161	ZyeloLight SPEC BIRC3/MAL T1 Dual Color Fusion Probe	11q22.1-q22.2 BIRC3 GF /1q21.31-q21.32 MAL T1 OF: DCD	Z-2146-50	0.05 ml	1	5%	250.00	262.50	262.5	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2146</a>
162	ZyeloLight SPEC RET Dual Color Break Apart Probe	10q11.21 RET DC BA: OF/GF	Z-2148-200	0.2 ml	1	5%	925.00	971.25	971.25	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2148</a>
163	ZyeloLight SPEC RET Dual Color Break Apart Probe	10q11.21 RET DC BA: OF/GF	Z-2148-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2148</a>
164	ZyeloLight SPEC USP6 Dual Color Break Apart Probe	17p13.2 USP6 DC BA: OF/GF	Z-2151-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2151</a>
165	ZyeloLight SPEC RREB1/MYB/CEN 6 Triple Color Probe	RREB1/MYB/CEN 6 Triple Color Probe	Z-2152-200	0.2 ml	1	5%	1331.00	1397.55	1397.55	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2152</a>
166	ZyeloLight SPEC RREB1/MYB/CEN 6 Triple Color Probe	RREB1/MYB/CEN 6 Triple Color Probe	Z-2152-50	0.05 ml	1	5%	400.00	420.00	420	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2152</a>
167	ZyeloLight SPEC TP53/CEN 17 Dual Color Probe	17p13.1 TP53 OF/CEN 17 (D17Z1) GF: DC	Z-2153-200	0.2 ml	1	5%	850.00	892.50	892.5	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2153</a>
168	ZyeloLight SPEC TP53/CEN 17 Dual Color Probe	17p13.1 TP53 OF/CEN 17 (D17Z1) GF: DC	Z-2153-50	0.05 ml	1	5%	250.00	262.50	262.5	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2153</a>
169	ZyeloLight SPEC ETV6/RUNX1 Dual Color Fusion Probe	ETV6 OF /RUNX1 GF: DC DF	Z-2157-200	0.2 ml	1	5%	780.00	819.00	819	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2157</a>
170	ZyeloLight SPEC ETV6/RUNX1 Dual Color Fusion Probe	ETV6 OF /RUNX1 GF: DC DF	Z-2157-50	0.05 ml	1	5%	235.00	246.75	246.75	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2157</a>
171	ZyeloLight SPEC TP53/ATM Dual Color Probe	17p13.1 TP53 OF/11q22.3 ATM GF: DC	Z-2159-200	0.2 ml	1	5%	820.00	861.00	861	<a href="https://www.zylovision.com/pro">https://www.zylovision.com/pro</a> <a href="https://www.zylovision.com/pro">ducts/zyolight/z-2159</a>

172	ZyeloLight SPEC TP59/ATM Dual Color Probe	17p13.1 TP59 OF/11q22.3 ATM GF, DC	Z-2159-50	0.05 ml	1	5%	250.00	262.50	262.5	<a href="https://www.zylovision.com/products/zyolight/z-2159">https://www.zylovision.com/products/zyolight/z-2159</a>
173	ZyeloLight SPEC D13S319/13q34/CEN 12 Triple Color Probe	13q14.2 D13S319 OF/13q34 BF /CEN 12 (D1223) GF, TC	Z-2160-200	0.2 ml	1	5%	1200.00	1260.00	1260	<a href="https://www.zylovision.com/products/zyolight/z-2160">https://www.zylovision.com/products/zyolight/z-2160</a>
174	ZyeloLight SPEC D13S319/13q34/CEN 12 Triple Color Probe	13q14.2 D13S319 OF/13q34 BF /CEN 12 (D1223) GF, TC	Z-2160-50	0.05 ml	1	5%	365.00	372.75	372.75	<a href="https://www.zylovision.com/products/zyolight/z-2160">https://www.zylovision.com/products/zyolight/z-2160</a>
175	ZyeloLight SPEC ALK2q11 Dual Color Probe	2p23.1-p23.2 ALK GF/2q11 OF, DC	Z-2161-200	0.2 ml	1	5%	950.00	997.50	997.5	<a href="https://www.zylovision.com/products/zyolight/z-2161">https://www.zylovision.com/products/zyolight/z-2161</a>
176	ZyeloLight SPEC ROS1/CEN 6 Dual Color Probe	6q22.1 ROS1 GF/CEN 6 (D6Z1) OF, DC	Z-2162-200	0.2 ml	1	5%	950.00	997.50	997.5	<a href="https://www.zylovision.com/products/zyolight/z-2162">https://www.zylovision.com/products/zyolight/z-2162</a>
177	ZyeloLight SPEC 18/CEN XY Triple Color Probe	DXZ1/DYZ3 18q21.31-q21.32/XP11.1-q11.1/YP11.1-q11.1	Z-2163-200	0.2 ml	1	5%	800.00	840.00	840	<a href="https://www.zylovision.com/products/zyolight/z-2163">https://www.zylovision.com/products/zyolight/z-2163</a>
178	ZyeloLight SPEC 13/21 Dual Color Probe	13q12.11/21q22.13-q22.2	Z-2164-200	0.2 ml	1	5%	450.00	472.50	472.5	<a href="https://www.zylovision.com/products/zyolight/z-2164">https://www.zylovision.com/products/zyolight/z-2164</a>
179	ZyeloLight SPEC RB1/13q12 Dual Color Probe	13q14.2 RB1 OF/13q12 GF 13q12.11, DC	Z-2165-200	0.2 ml	1	5%	870.00	913.50	913.5	<a href="https://www.zylovision.com/products/zyolight/z-2165">https://www.zylovision.com/products/zyolight/z-2165</a>
180	ZyeloLight SPEC RB1/13q12 Dual Color Probe	13q14.2 RB1 OF/13q12 GF 13q12.11, DC	Z-2165-50	0.05 ml	1	5%	260.00	273.00	273	<a href="https://www.zylovision.com/products/zyolight/z-2165">https://www.zylovision.com/products/zyolight/z-2165</a>
181	FlexISH ERBB2/CEN 17 Dual Color Probe	ERBB2/CEN 17 DC, GF/OF	Z-2166-200	0.2 ml	1	5%	1100.00	1155.00	1155	<a href="https://www.zylovision.com/products/zyolight/z-2166">https://www.zylovision.com/products/zyolight/z-2166</a>
182	FlexISH ERBB2/CEN 17 Dual Color Probe	ERBB2/CEN 17 DC, GF/OF	Z-2166-50	0.05 ml	1	5%	345.00	362.25	362.25	<a href="https://www.zylovision.com/products/zyolight/z-2166">https://www.zylovision.com/products/zyolight/z-2166</a>
183	ZyeloLight SPEC NTRK1 Dual Color Break Apart Probe	1q22-q23 NTRK1 DC BA: GF/OF	Z-2167-200	0.2 ml	1	5%	900.00	945.00	945	<a href="https://www.zylovision.com/products/zyolight/z-2167">https://www.zylovision.com/products/zyolight/z-2167</a>
184	ZyeloLight SPEC NTRK1 Dual Color Break Apart Probe	1q22-q23 NTRK1 DC BA: GF/OF	Z-2167-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/products/zyolight/z-2167">https://www.zylovision.com/products/zyolight/z-2167</a>
185	ZyeloLight SPEC FGFR1 Dual Color Break Apart Probe	8p11.23-p11.22 FGFR1 DC BA: OF/GF	Z-2168-200	0.2 ml	1	5%	870.00	913.50	913.5	<a href="https://www.zylovision.com/products/zyolight/z-2168">https://www.zylovision.com/products/zyolight/z-2168</a>
186	ZyeloLight SPEC FGFR1 Dual Color Break Apart Probe	8p11.23-p11.22 FGFR1 DC BA: OF/GF	Z-2168-50	0.05 ml	1	5%	270.00	283.50	283.5	<a href="https://www.zylovision.com/products/zyolight/z-2168">https://www.zylovision.com/products/zyolight/z-2168</a>
187	ZyeloLight SPEC FGFR2 Dual Color Break Apart Probe	10q26.12-q26.13 FGFR2 DC BA: OF/GF	Z-2169-200	0.2 ml	1	5%	900.00	945.00	945	<a href="https://www.zylovision.com/products/zyolight/z-2169">https://www.zylovision.com/products/zyolight/z-2169</a>
188	ZyeloLight SPEC FGFR2 Dual Color Break Apart Probe	4p16.3 FGFR3 DC BA: OF/GF	Z-2170-200	0.2 ml	1	5%	870.00	913.50	913.5	<a href="https://www.zylovision.com/products/zyolight/z-2170">https://www.zylovision.com/products/zyolight/z-2170</a>
189	ZyeloLight SPEC FGFR3 Dual Color Break Apart Probe	4p16.3 FGFR3 DC BA: OF/GF	Z-2170-50	0.05 ml	1	5%	270.00	283.50	283.5	<a href="https://www.zylovision.com/products/zyolight/z-2170">https://www.zylovision.com/products/zyolight/z-2170</a>
190	ZyeloLight SPEC BCL2L1/CEN 20 Dual Color Probe	BCL2L1/CEN 20 DC	Z-2171-200	0.2 ml	1	5%	900.00	945.00	945	<a href="https://www.zylovision.com/products/zyolight/z-2171">https://www.zylovision.com/products/zyolight/z-2171</a>
191	ZyeloLight SPEC MCL1/1p12 Dual Color Probe	MCL1/1p12 DC	Z-2173-200	0.2 ml	1	5%	900.00	945.00	945	<a href="https://www.zylovision.com/products/zyolight/z-2173">https://www.zylovision.com/products/zyolight/z-2173</a>
192	ZyeloLight SPEC BCL2/CEN 18 Dual Color Probe	18q21.33 BCL2 GF/CEN 18 (D18Z1) OF, DC	Z-2174-50	0.05 ml	1	5%	270.00	283.50	283.5	<a href="https://www.zylovision.com/products/zyolight/z-2174">https://www.zylovision.com/products/zyolight/z-2174</a>
193	ZyeloLight SPEC YWHA6 Dual Color Break Apart Probe	YWHA6 DC BA	Z-2175-50	0.05 ml	8	5%	284.00	298.20	2385.6	<a href="https://www.zylovision.com/products/zyolight/z-2175">https://www.zylovision.com/products/zyolight/z-2175</a>
194	ZyeloLight SPEC ETV6 Dual Color Break Apart Probe	12p13.2 ETV6 DC BA: GF/OF	Z-2176-200	0.2 ml	1	5%	810.00	850.50	850.5	<a href="https://www.zylovision.com/products/zyolight/z-2176">https://www.zylovision.com/products/zyolight/z-2176</a>
195	ZyeloLight SPEC ETV6 Dual Color Break Apart Probe	12p13.2 ETV6 DC BA: GF/OF	Z-2176-50	0.05 ml	1	5%	260.00	273.00	273	<a href="https://www.zylovision.com/products/zyolight/z-2176">https://www.zylovision.com/products/zyolight/z-2176</a>
196	ZyeloLight SPEC BCL6 Dual Color Break Apart Probe	3q27.3-q28 BCL6 DC BA: GF/OF	Z-2177-200	0.2 ml	1	5%	1090.00	1144.50	1144.5	<a href="https://www.zylovision.com/products/zyolight/z-2177">https://www.zylovision.com/products/zyolight/z-2177</a>
197	ZyeloLight SPEC BCL6 Dual Color Break Apart Probe	3q27.3-q28 BCL6 DC BA: GF/OF	Z-2177-50	0.05 ml	1	5%	350.00	367.50	367.5	<a href="https://www.zylovision.com/products/zyolight/z-2177">https://www.zylovision.com/products/zyolight/z-2177</a>
198	ZyeloLight SPEC SMARCB1/22q12 Dual Color Probe	SMARCB1/22q12 DC	Z-2178-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zylovision.com/products/zyolight/z-2178">https://www.zylovision.com/products/zyolight/z-2178</a>

199	ZyLight SPEC CD274, PDCC1LG2/CEN 9 Dual Color Probe	CD274, PDCC1LG2/CEN 9 DC	Z-2179-200	0.2 ml	1	5%	900,00	945,00	945	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2179">https://www.zyovision.com/p/roducts/zytolight/z-2179</a>
200	ZyLight SPEC CD274, PDCC1LG2/CEN 9 Dual Color Probe	CD274, PDCC1LG2/CEN 9 DC	Z-2179-50	0.05 ml	1	5%	280,00	294,00	294	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2179">https://www.zyovision.com/p/roducts/zytolight/z-2179</a>
201	ZyLight SPEC 21/CEN XYq12 Triple Color Probe	DXZ1/III DY21 21q22.13-q22.2/Xp11.1-q11.1/Yq12	Z-2180-200	0.2 ml	1	5%	800,00	840,00	840	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2180">https://www.zyovision.com/p/roducts/zytolight/z-2180</a>
202	ZyLight SPEC NRG1 Dual Color Break Apart Probe	NRG1 DC BA	Z-2181-200	0.2 ml	1	5%	915,00	960,75	960,75	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2181">https://www.zyovision.com/p/roducts/zytolight/z-2181</a>
203	FlexISH-Tissue Implementation Kit	Heat Pretreatment Solution Citric, 500 ml; Pepsin Solution, 4 ml; 5x FlexISH Wash Buffer, 500 ml; DAPI/DuratTect-Solution, 0.8 ml	Z-2182-20	20 testu	1	5%	240,00	252,00	252	<a href="https://www.zyovision.com/p/roducts/flexish">https://www.zyovision.com/p/roducts/flexish</a>
204	FlexISH-Tissue Implementation Kit	Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; 5x FlexISH Wash Buffer, 150 ml; DAPI/DuratTect-Solution, 0.2 ml	Z-2182-5	5 testai	1	5%	100,00	105,00	105	<a href="https://www.zyovision.com/p/roducts/flexish">https://www.zyovision.com/p/roducts/flexish</a>
205	ZyLight SPEC EWSR1/FL11 TrnCheck Probe	EWSR1/FL11 TrnCheck Probe	Z-2183-50	0.05 ml	1	5%	320,00	336,00	336	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2183">https://www.zyovision.com/p/roducts/zytolight/z-2183</a>
206	ZyLight SPEC SS18/SSX1 TrnCheck Probe	SS18/SSX1 TrnCheck Probe	Z-2184-50	0.05 ml	1	5%	320,00	336,00	336	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2184">https://www.zyovision.com/p/roducts/zytolight/z-2184</a>
207	ZyLight SPEC FOXO1/PAX3 TrnCheck Probe	FOXO1/PAX3 TrnCheck Probe	Z-2185-50	0.05 ml	1	5%	370,00	388,50	388,5	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2185">https://www.zyovision.com/p/roducts/zytolight/z-2185</a>
208	ZyLight SPEC BRAF Dual Color Break Apart Probe	FGFR3 BAF/GFR2/FGFR1 Quadruple Color Probe	Z-2189-200	0.2 ml	1	5%	920,00	966,00	966	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2189">https://www.zyovision.com/p/roducts/zytolight/z-2189</a>
209	ZyLight SPEC ERBB2/17S122 Dual Color Probe	ERBB2/17S122 DC	Z-2190-200	0.2 ml	1	5%	1000,00	1050,00	1050	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2190">https://www.zyovision.com/p/roducts/zytolight/z-2190</a>
210	ZyLight SPEC ERBB2/17S122 Dual Color Probe	ERBB2/17S122 DC	Z-2190-50	0.05 ml	1	5%	305,00	320,25	320,25	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2190">https://www.zyovision.com/p/roducts/zytolight/z-2190</a>
211	ZyLight SPEC BRAF/CEN 7 Dual Color Probe	BRAF/CEN 7 DC	Z-2191-200	0.2 ml	1	5%	1000,00	1050,00	1050	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2191">https://www.zyovision.com/p/roducts/zytolight/z-2191</a>
212	ZyLight SPEC BCL2 Dual Color Break Apart Probe	BCL2DC BA	Z-2192-200	0.2 ml	1	5%	1300,00	1365,00	1365	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2192">https://www.zyovision.com/p/roducts/zytolight/z-2192</a>
213	ZyLight SPEC BCL2 Dual Color Break Apart Probe	BCL2DC BA	Z-2192-50	0.05 ml	1	5%	400,00	420,00	420	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2192">https://www.zyovision.com/p/roducts/zytolight/z-2192</a>
214	ZyLight SPEC KMT2A Dual Color Break Apart Probe	KMT2A DC BA	Z-2193-200	0.2 ml	1	5%	790,00	829,50	829,5	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2193">https://www.zyovision.com/p/roducts/zytolight/z-2193</a>
215	ZyLight SPEC KMT2A Dual Color Break Apart Probe	KMT2A DC BA	Z-2193-50	0.05 ml	1	5%	250,00	262,50	262,5	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2193">https://www.zyovision.com/p/roducts/zytolight/z-2193</a>
216	ZyLight SPEC NRG1/CD74 TrnCheck Probe	NRG1/CD74 TrnCheck Probe	Z-2194-200	0.2 ml	1	5%	1300,00	1365,00	1365	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2194">https://www.zyovision.com/p/roducts/zytolight/z-2194</a>
217	ZyLight SPEC VEGFACEN 6 Dual Color Probe	VEGFACEN 6 DC	Z-2195-200	0.2 ml	1	5%	950,00	997,50	997,5	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2195">https://www.zyovision.com/p/roducts/zytolight/z-2195</a>
218	ZyLight SPEC MALT1 Dual Color Break Apart Probe	MALT1 DC BA	Z-2196-200	0.2 ml	1	5%	770,00	808,50	808,5	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2196">https://www.zyovision.com/p/roducts/zytolight/z-2196</a>
219	ZyLight SPEC MALT1 Dual Color Break Apart Probe	MALT1 DC BA	Z-2196-50	0.05 ml	1	5%	250,00	262,50	262,5	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2196">https://www.zyovision.com/p/roducts/zytolight/z-2196</a>
220	ZyLight SPEC PDGFRB Dual Color Break Apart Probe	PDGFRB DC BA	Z-2197-50	0.05 ml	1	5%	250,00	262,50	262,5	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2197">https://www.zyovision.com/p/roducts/zytolight/z-2197</a>
221	ZyLight SPEC TP53/17q22 Dual Color Probe	TP53/17q22 DC	Z-2198-50	0.05 ml	1	5%	260,00	273,00	273	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2198">https://www.zyovision.com/p/roducts/zytolight/z-2198</a>
222	ZyLight SPEC ABL1 Dual Color Break Apart Probe	ABL1 DC BA	Z-2199-50	0.05 ml	1	5%	260,00	273,00	273	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2199">https://www.zyovision.com/p/roducts/zytolight/z-2199</a>
223	ZyLight SPEC ABL2 Dual Color Break Apart Probe	ABL2 DC Bap	Z-2200-50	0.05 ml	1	5%	260,00	273,00	273	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2200">https://www.zyovision.com/p/roducts/zytolight/z-2200</a>
224	ZyLight SPEC CRLF2 Dual Color Break Apart Probe	CRLF2 DC BA	Z-2201-50	0.05 ml	1	5%	260,00	273,00	273	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2201">https://www.zyovision.com/p/roducts/zytolight/z-2201</a>
225	ZyLight SPEC CSF1R Dual Color Break Apart Probe	CSF1R DC BA	Z-2202-50	0.05 ml	1	5%	260,00	273,00	273	<a href="https://www.zyovision.com/p/roducts/zytolight/z-2202">https://www.zyovision.com/p/roducts/zytolight/z-2202</a>

226	FlexiSH ALK/ROS1 Distinguish Probe	ALK/ROS1 Distinguish Probe	Z-2203-200	0.2 ml	1	5%	1300.00	1365.00	1365	<a href="https://www.zyovision.com/products/flexish/z-2203">https://www.zyovision.com/products/flexish/z-2203</a>
227	FlexiSH ALK/ROS1 Distinguish Probe	ALK/ROS1 Distinguish Probe	Z-2203-50	0.05 ml	1	5%	390.00	409.50	409.5	<a href="https://www.zyovision.com/products/flexish/z-2203">https://www.zyovision.com/products/flexish/z-2203</a>
228	ZyLight SPEC NTRK2 Dual Color Break Apart Probe	NTRK2 DC BA	Z-2205-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zyovision.com/products/zyolight/z-2205">https://www.zyovision.com/products/zyolight/z-2205</a>
229	ZyLight SPEC NTRK2 Dual Color Break Apart Probe	NTRK2 DC BA	Z-2205-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zyovision.com/products/zyolight/z-2205">https://www.zyovision.com/products/zyolight/z-2205</a>
230	ZyLight SPEC NTRK3 Dual Color Break Apart Probe	NTRK3 DC BA	Z-2206-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zyovision.com/products/zyolight/z-2206">https://www.zyovision.com/products/zyolight/z-2206</a>
231	ZyLight SPEC NTRK3 Dual Color Break Apart Probe	NTRK3 DC BA	Z-2206-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zyovision.com/products/zyolight/z-2206">https://www.zyovision.com/products/zyolight/z-2206</a>
232	ZyLight SPEC CBF3 Dual Color Break Apart Probe	CBFB DC BA	Z-2207-50	0.05 ml	1	5%	280.00	294.00	294	<a href="https://www.zyovision.com/products/zyolight/z-2207">https://www.zyovision.com/products/zyolight/z-2207</a>
233	ZyLight SPEC NUTM1 Dual Color Break Apart Probe	NUTM1 DC BA	Z-2208-200	0.2 ml	1	5%	920.00	966.00	966	<a href="https://www.zyovision.com/products/zyolight/z-2208">https://www.zyovision.com/products/zyolight/z-2208</a>
234	ZyLight SPEC PDGFRA/FIP1L1 TriCheck Probe	PDGFRA/FIP1L1 TriCheck Probe	Z-2209-50	0.05 ml	1	5%	330.00	346.50	346.5	<a href="https://www.zyovision.com/products/zyolight/z-2209">https://www.zyovision.com/products/zyolight/z-2209</a>
235	ZyLight SPEC IRF4/DUSP22 Dual Color Break Apart Probe	IRF4/DUSP22 DC BA	Z-2210-50	0.05 ml	8	5%	347.00	364.35	2914.8	<a href="https://www.zyovision.com/products/zyolight/z-2210">https://www.zyovision.com/products/zyolight/z-2210</a>
236	ZyLight SPEC ESR1/D5S23_D5S721 Dual Color Probe	ESR1/D5S23_D5S721 DC	Z-2211-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/products/zyolight/z-2211">https://www.zyovision.com/products/zyolight/z-2211</a>
237	ZyLight SPEC WWTR1 Dual Color Break Apart Probe	WWTR1 DC BA	Z-2212-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/products/zyolight/z-2212">https://www.zyovision.com/products/zyolight/z-2212</a>
238	ZyLight SPEC PTPRT/20q11 Dual Color Probe	PTPRT/20q11 DC	Z-2213-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/products/zyolight/z-2213">https://www.zyovision.com/products/zyolight/z-2213</a>
239	ZyLight SPEC CUX1/EZH2/CEN 7 Triple Color Probe	CUX1/EZH2/CEN 7 TC	Z-2214-50	0.05 ml	1	5%	320.00	336.00	336	<a href="https://www.zyovision.com/products/zyolight/z-2214">https://www.zyovision.com/products/zyolight/z-2214</a>
240	ZyLight SPEC PHF1 Dual Color Break Apart Probe	PHF1 DC BA	Z-2215-50	0.05 ml	1	5%	260.00	262.50	262.5	<a href="https://www.zyovision.com/products/zyolight/z-2215">https://www.zyovision.com/products/zyolight/z-2215</a>
241	ZyLight SPEC 11q gainless Triple Color Probe	11q gainless TC	Z-2216-50	0.05 ml	2	5%	347.00	364.35	728.7	<a href="https://www.zyovision.com/products/zyolight/z-2216">https://www.zyovision.com/products/zyolight/z-2216</a>
242	ZyLight SPEC NUP214 Dual Color Break Apart Probe	NUP214 DC BA	Z-2266-50	0.05 ml	1	5%	220.00	231.00	231	<a href="https://www.zyovision.com/products/zyolight/z-2266">https://www.zyovision.com/products/zyolight/z-2266</a>
243	ZyLight SPEC NUP98 Dual Color Break Apart Probe	NUP98 DC BA	Z-2266-50	0.05 ml	1	5%	220.00	231.00	231	<a href="https://www.zyovision.com/products/zyolight/z-2266">https://www.zyovision.com/products/zyolight/z-2266</a>
244	ZyLight SPEC CREBBP Dual Color Break Apart Probe	CREBBP DC BA	Z-2267-50	0.05 ml	1	5%	220.00	231.00	231	<a href="https://www.zyovision.com/products/zyolight/z-2267">https://www.zyovision.com/products/zyolight/z-2267</a>
245	ZyLight SPEC CSF1R/D5S23_D5S721 Dual Color Probe	CSF1R/D5S23_D5S721 DC	Z-2268-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/products/zyolight/z-2268">https://www.zyovision.com/products/zyolight/z-2268</a>
246	FlexiSH RET/KIF5B TriCheck Probe	RET 10q11.21;KIF5B 10q11.22; TC; Green/Orange/Blue	Z-2269-200	0.2 ml	1	5%	1100.00	1155.00	1155	<a href="https://www.zyovision.com/products/flexish/z-2269">https://www.zyovision.com/products/flexish/z-2269</a>
247	FlexiSH RET/KIF5B TriCheck Probe	RET 10q11.21;KIF5B 10q11.22; TC; Green/Orange/Blue	Z-2269-50	0.05 ml	1	5%	330.00	367.50	367.5	<a href="https://www.zyovision.com/products/flexish/z-2269">https://www.zyovision.com/products/flexish/z-2269</a>
248	ZyLight SPEC MAF/IGH Dual Color Dual Fusion Probe	MAF/IGH DC DF	Z-2270-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/products/zyolight/z-2270">https://www.zyovision.com/products/zyolight/z-2270</a>
249	ZyLight SPEC MAFB/IGH Dual Color Dual Fusion Probe	MAFB/IGH DC DF	Z-2271-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/products/zyolight/z-2271">https://www.zyovision.com/products/zyolight/z-2271</a>
250	ZyLight Glioma 1p/19q Probe Set	Glioma 1p/19q	Z-2272-20	2x 0.2 ml	5	5%	1575.00	1653.75	8268.75	<a href="https://www.zyovision.com/products/zyolight/z-2272">https://www.zyovision.com/products/zyolight/z-2272</a>
251	ZyLight SPEC TERT Dual Color Break Apart Probe	TERT DC BA	Z-2273-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/products/zyolight/z-2273">https://www.zyovision.com/products/zyolight/z-2273</a>
252	ZyLight SPEC C19MC/19p13 Dual Color Probe	C19MC/19p13 DC	Z-2274-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/products/zyolight/z-2274">https://www.zyovision.com/products/zyolight/z-2274</a>

263	ZyloLight SPEC ZNF384 Dual Color Break Apart Probe	ZNF384 DC BA	Z-2275-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2275">https://www.zyovision.com/price/ducts/zyolight/z-2275</a>
264	ZyloLight SPEC CKS1B/CDKN2C Dual Color Probe	CKS1B/CDKN2C DC	Z-2276-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2276">https://www.zyovision.com/price/ducts/zyolight/z-2276</a>
265	ZyloLight SPEC MEF2D/BCL9 TriCheck Probe	MEF2D/BCL9 TriCheck Probe	Z-2277-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2277">https://www.zyovision.com/price/ducts/zyolight/z-2277</a>
266	ZyloLight SPEC RICTOR/5q31.1 Dual Color Probe	RICTOR/5q31.1 DC, Green/Orange	Z-2278-200	0.2 ml	1	5%	950.00	997.50	997.5	<a href="https://www.zyovision.com/price/zyolight/price-index">https://www.zyovision.com/price/zyolight/price-index</a>
267	ZyloLight Aneuploidy Panel 18/XY and 13/21	18/XY and 13/21	Z-2279-20	2x 0.2 ml	1	5%	1100.00	1155.00	1155	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2279">https://www.zyovision.com/price/ducts/zyolight/z-2279</a>
268	ZyloLight SPEC D13S319/13q34 Dual Color Probe	D13S319/13q34 DC	Z-2280-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2280">https://www.zyovision.com/price/ducts/zyolight/z-2280</a>
269	ZyloLight SPEC MYB/CEN 6 Dual Color Probe	MYB/CEN 6 DC	Z-2281-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2281">https://www.zyovision.com/price/ducts/zyolight/z-2281</a>
270	ZyloLight SPEC FGFR3/IGH Dual Color Dual Fusion Probe	FGFR3/IGH DC	Z-2282-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2282">https://www.zyovision.com/price/ducts/zyolight/z-2282</a>
271	FlexISH BCL2/BCL6 DistinguISH Probe	BCL2/BCL6 DistinguISH Probe	Z-2283-200	0.2 ml	1	5%	1400.00	1470.00	1470	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2283">https://www.zyovision.com/price/ducts/zyolight/z-2283</a>
272	FlexISH BCL2/BCL6 DistinguISH Probe	BCL2/BCL6 DistinguISH Prob	Z-2283-50	0.05 ml	1	5%	400.00	420.00	420	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2283">https://www.zyovision.com/price/ducts/zyolight/z-2283</a>
273	ZyloLight SPEC TERC/CEN 3 Dual Color Probe	TERC/CEN 3 DC	Z-2284-200	0.2 ml	1	5%	950.00	997.50	997.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2284">https://www.zyovision.com/price/ducts/zyolight/z-2284</a>
274	ZyloLight SPEC CIC Dual Color Break Apart Probe	CIC DC BA	Z-2285-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2285">https://www.zyovision.com/price/ducts/zyolight/z-2285</a>
275	ZyloLight SPEC IGL Dual Color Break Apart Probe	IGL DC BA	Z-2286-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2286">https://www.zyovision.com/price/ducts/zyolight/z-2286</a>
276	ZyloLight SPEC GATA2/MECOM1 Dual Color Dual Fusion Probe	GATA2/MECOM DC DF, Orange/Green	Z-2287-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2287">https://www.zyovision.com/price/ducts/zyolight/z-2287</a>
277	ZyloLight SPEC IGH Dual Color Break Apart Probe	IGH DC BA, Green/Orange	Z-2288-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2288">https://www.zyovision.com/price/ducts/zyolight/z-2288</a>
278	ZyloLight SPEC D1George Triple Color Probe	D1George TC, Blue/Orange/Green	Z-2289-50	0.05 ml	1	5%	200.00	210.00	210	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2289">https://www.zyovision.com/price/ducts/zyolight/z-2289</a>
279	ZyloLight SPEC SPI1 Dual Color Break Apart Probe	SPI1 DC BA, Green/Orange	Z-2291-50	0.05 ml	1	5%	230.00	241.50	241.5	<a href="https://www.zyovision.com/price/ducts/zyolight/z-2291">https://www.zyovision.com/price/ducts/zyolight/z-2291</a>
280	ZyloMation ERBB2/CEN 17 Dual Color FISH Probe	ERBB2/CEN 17 DC, Green/Orange	Z-2292-5 1ML	5.1 ml	1	5%	1460.00	1533.00	1533	<a href="https://www.zyovision.com/products/zyomation/z-2292">https://www.zyovision.com/products/zyomation/z-2292</a>
281	ZyloMation ROS1 Dual Color Break Apart FISH Probe	ROS1 DC BA, Green/Orange	Z-2298-5 1ML	5.1 ml	1	5%	1460.00	1533.00	1533	<a href="https://www.zyovision.com/products/zyomation/z-2298">https://www.zyovision.com/products/zyomation/z-2298</a>
282	ZyloFast PLUS CISH Implementation Kit	AP-NBT/BCIP	T-1061-40	mkl (40 test)	1	5%	430.00	451.50	451.5	<a href="https://www.zyovision.com/products/zyofast/z-1061">https://www.zyovision.com/products/zyofast/z-1061</a>
283	ZyloFast CMV Probe	Digoxigenin	T-1113-400	400 mkl	1	5%	115.00	120.75	120.75	<a href="https://www.zyovision.com/products/zyofast/z-1113">https://www.zyovision.com/products/zyofast/z-1113</a>
284	ZyloFast EBV Probe	Digoxigenin	T-1114-400	400 mkl	1	5%	115.00	120.75	120.75	<a href="https://www.zyovision.com/products/zyofast/z-1114">https://www.zyovision.com/products/zyofast/z-1114</a>

275	ZytoFast human Ig-kappa Probe	Digoxigenin	T-1115-400	400 mkl	1	5%	115,00	120,75	120,75	<a href="https://www.zytovision.com/products/z-yofastplus-t-1115">https://www.zytovision.com/products/z-yofastplus-t-1115</a>
276	ZytoFast human Ig-lambda Probe	Digoxigenin	T-1116-400	400 mkl	1	5%	115,00	120,75	120,75	<a href="https://www.zytovision.com/products/z-yofastplus-t-1116">https://www.zytovision.com/products/z-yofastplus-t-1116</a>

**Pardavėjas**

UAB „Medfarmos laboratorijos“

P. Smuglevičiaus g. 1, Vilnius LT-08311

Įmonės kodas 124484379

Direktorė

Živilė Sukackienė

A.V.




**Pirkėjas**  
 VšĮ Vilniaus universiteto ligoninė Santaros klinikos  
 Santariškių g. 2, LT-08661 Vilnius  
 Įmonės kodas 124364561  
 PVM mok. kodas LT243645610  
 A. s. LT71 7300 0100 0249 2260  
 AB „Swedbank“ b. k. 73000  
 Tel. (8 5) 236 5000, faks. (8 5) 236 5111

Generalinis direktorius  
 Feliksas Jankevičius  
 A.V.



