




TKBAG[®]

SPECIMEN RETRIEVAL BAG



Safety

Extra-durable bag prevents leaks and tears. Soft ring design reduces tissue damage risks.

Convenience: Easy 4 U.

Easy in, easy out, easy open, easy close; Saves time! Multiple sizes including extra-large

Pouch with string is optional

Quality

Strict, extensive raw material selection and 31 quality inspection steps on every TKBAG.

The TKBAG/ Retriever Bag is used in all types of laparoscopic surgery for removing excised tissues/specimen through a small surgical incision. It can be introduced through 10mm or larger size trocar canula.

It comes in a variety of sizes, (as shown in the table below) and is suitable for a wide range of procedures including, Gallbladders, Appendices, Ectopic pregnancies, Lymph nodes, Ovaries, Lung wedges, Sections of bowel, Splenectomies, Nephrectomies and more.

Product Code	Bag Diam (mm)	Bag length (mm)	Volume (ml)	Qty/Box
HSD-60	60	120	75	10
HSD-80	80	175	200	10
HSD-80/C	80	175	200	10
HSD-100	100	170	350	10
HSD-130	130	175	700	10



CE 0197
ISO9001
ISO13485